CVS/pharmacy
by **microlife**

DIGITAL ADVANCED AUTOMATIC

Blood Pressure Monitor

Model #BP3AC1-1PC
ID #3AC1-AP/271244

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American Heart
Association®

Learn and Live.

MicroLife Corporation is a proud national
sponsor of the American Heart Association's
High Blood Pressure Heart Profilers™

PLEASE NOTE:
THIS MEDICAL INSTRUMENT MUST BE USED ACCORDING
TO INSTRUCTIONS TO ENSURE ACCURATE READINGS.

Automatic Blood Pressure Monitor

Instruction Manual

Table of contents

1. Introduction

- 1.1. Features of your blood pressure monitor
- 1.2. Important information about self measurement

2. Important information on blood pressure and its measurement

- 2.1. How does high/low blood pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done if regular high/low values are obtained?
- 2.4. Irregular Heartbeat Detector
- 2.5. MAM technology

3. Components of your blood pressure monitor

4. Operation of your blood pressure monitor

- 4.1. Inserting the batteries
- 4.2. Using an A/C power adapter
- 4.3. Tube connection

5. Measurement Procedure

- 5.1. Before measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Setting the time and date
- 5.5. Reading the set date
- 5.6. Selection of Microlife Average Mode
- 5.7. Selecting the User
- 5.8. Measuring procedure
- 5.9. Memory - Displaying the last 99 measurements
- 5.10. Discontinuing a measurement

6. Error Messages/Troubleshooting

7. Care and maintenance

8. Warranty

9. Certifications

10. Technical specifications

11. How to contact us

1. Introduction

1.1. Features of your Digital Advanced Automatic blood pressure monitor,

Your CVS/Microlife blood pressure monitor with MAM (microlife Averaging Mode) Technology and Irregular Heartbeat Detection is a fully automatic digital blood pressure measuring device for use by adults on the upper arm at home or in your doctor's/nurse's office. It enables very fast and reliable measurement of the systolic and diastolic blood pressure as well as the pulse by way of the oscillometric method. This device offers clinically proven accuracy and has been designed to be user friendly.

Before using, please read this instruction manual carefully and then keep it in a safe place. Please contact your doctor for further questions on the subject of blood pressure and its measurement.

Attention!

1.2. Important information about self-measurement

- **Self-measurement means Control**, not diagnosis or treatment. Unusual values must always be discussed with your doctor. **Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.**
- The pulse display is **not** suitable for checking the frequency of heart pacemakers!
- In cases of irregular heartbeat, measurements made with this instrument should only be evaluated after consultation with your doctor.

Electromagnetic interference:

The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring accuracy.

Blood Pressure Analyzer Software (Sold Separately as an Optional Accessory):

This unit can be used in connection with your PC using the Microlife Blood Pressure Analyzer software and a USB cable. Your PC will allow a capacity of monitoring 80 patients, each with 1000 data (note: overuse will lower system efficiency). For information about the software and how to purchase, call 1-866-464-6184 or visit www.microlifeusa.com.

System Requirements for Blood Pressure Analyzer Software:

- Windows 98 SE, 2000, XP
- CD-Rom Drive
- Minimum 256MB RAM
- 500MB Available Hard Disk Space
- USB Port version 1.0 or higher
- Microlife Blood Pressure Monitor with USB Port

2. Important information on blood pressure and its measurement

2.1. How does high/low blood pressure arise?

Your level of blood pressure is determined in the circulatory center of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is effected by fine muscles in the blood vessel walls.

Your level of arterial blood pressure changes periodically during heart activity: During the “blood ejection” (Systole) the value is highest (systolic blood pressure value). At the end of the heart’s “rest period” (Diastole) pressure is lowest (diastolic blood pressure value).

Blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2. Which values are normal?

Blood pressure is too high if your diastolic pressure is above 90 mmHg and/or your systolic blood pressure is over 160 mmHg, **while at rest**. In this case, please consult your physician immediately. Long-term values at this level endanger your health due to continual damage to the blood vessels in you body.

If your systolic blood pressure values are between 140 mmHg and 160mmHg and/or the diastolic blood pressure values between 90 mmHg and 95mmHg, consult your physician. Regular self-checks are necessary.

If you have blood pressure values that are too low, (i.e., systolic values under 105mmHg and/or diastolic values under 60 mmHg), consult your physician.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your physician. **Never use the results of your measurements to independently alter the drug doses prescribed by your physician.**

Which values are normal?

The following standards for assessing high blood pressure (without regard to age) have been established by the National Institutes of Health JNC7, 2003.

Category	Systolic (mmHg)	Diastolic (mmHg)
Normal	<120	<80
Pre-Hypertension	120-139	80-89
Hypertension		
Stage 1	140-159	90-99
Stage 2	>160	>100

Further information

- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called “labile hypertension.” Consult your doctor.
- **Correctly measured diastolic blood pressure values above 120mmHg and require immediate medical treatment.**

2.3. What can be done if regular high or low values are obtained?

- a) Consult your doctor.
- b) Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result from arteriosclerosis. Furthermore, the heart will become structurally damaged with increased blood pressure values.
- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organ malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take to reduce and even prevent high blood pressure. These measures must be permanent lifestyle changes.

A) Eating habits

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

B) Previous illnesses

- Consistently follow all medical instructions for treating illness such as:
 - Diabetes (Diabetes mellitus or sugar diabetes)
 - Fat metabolism disorder
 - Gout

C) Habits

- Give up smoking completely.
- Drink only moderate amounts of alcohol.
- Restrict your caffeine consumption (coffee, tea, chocolate).

D) Your Physical condition

- **After a medical examination, and with your doctor's approval and direction; exercise.**
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your exercise routine. You must receive advise regarding the type and extent of exercise that is appropriate for you.

2.4. Irregular Heartbeat Detector

The appearance of this symbol  signifies that a certain irregular heartbeat was detected which could cause or indicate cardiovascular disease. This indicator is only a caution. The CVS/Microlife Fully Automatic Blood Pressure Monitor with Irregular Heartbeat Detector can detect blood pressure and pulse rate measurement with specifications even when the patient is having an irregular heartbeat. It is important that you be relaxed, remain still and do not talk during the measurements.

Note: We recommend contacting your physician if you see this “” icon frequently.


CAUTION: An irregular heartbeat can be a serious condition requiring medical attention. See your doctor if you suspect you have an irregular heartbeat.

2.5. MAM Technology

Microlife Averaging Mode (MAM) Technology is a new function for optimum reliability in self measurement of blood pressure. A single measurement is rarely your real result. An average of three readings is medically proven to provide the most accurate result. With the push of a single button, the Microlife Averaging Mode blood pressure monitor takes three consecutive measurements, and calculates the average.

3. Components of your blood pressure monitor


a) Measuring unit

b) Upper arm cuffs:

Type S101 for arm circumference 24–33.7 cm or
9.5" - 13.25" (included)

Type S102 for arm circumference 30.5 - 40.64 cm
12" - 16" (sold separately)

Please Note:

Arm circumference should be measured with a measuring tape in the middle of the relaxed upper arm. Do not force cuff connection into the opening. Make sure the cuff connection is not pushed into the AC adapter port. If the cuff is too small, call 1-866-464-6184 for further information. You may prefer to use one of our wrist type monitors.

4. Operation of your blood pressure monitor



4.1. Inserting the batteries

Insert batteries immediately after unpacking the unit. The battery compartment is located on the back of the device (see illustration).

a) Remove cover as illustrated

b) Insert the batteries (4 x size AA 1,5V),
observing the indicated polarity.


c) If a battery warning appears in the display,
the batteries are discharged and must be replaced.

4.2. Using an AC power adapter (special accessory)

It is possible to operate this blood pressure instrument with an AC adapter. (output 6 V DC / 600 mA with DIN plug). Make certain that you use an adapter which fulfills the legal requirements and electronic requirements in the U.S. (UL standard). Use of an AC adapter other than a Microlife or CVS adapter may invalidate your warranty.

- a) Push the plug into the socket on the back of the monitor.
- b) Plug the AC adapter into a 110 V power socket (U.S. or Canada).
- c) Test that power is available by pressing the On/Off button.


Note:

- No power is taken from the batteries while the AC adapter is connected to the instrument.
- If the power is interrupted during a measurement (e.g. by removal of the adapter from the wall socket), the instrument must be reset by removing the plug from the instrument.
- Please call 1-866-464-6184 if you have questions relating to the AC adapter.

4.3. Tube connection

Insert the cuff tube into the opening on the side of the monitor, as shown in the diagram.

5. Measurement Procedure

Please note: You should always be seated before and during measurement.

5.1. Before measurement:

- Avoid eating and smoking as well as all forms of exertion directly before measurement. These factors influence the measurement result. Find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Always perform measurements at the same time of day, since blood pressure changes during the course of the day.

5.2. Common sources of error:

Note: Comparable blood pressure measurements always require the same conditions! Conditions should always be quiet.

- All efforts by the user to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not flex any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower or higher than the heart, an erroneously high or low blood pressure will be measured! Each 15cm (6 inch) difference in height between your heart and the cuff results in a measurement error of 10 mmHg!
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is extremely important. Cuff size is dependent upon the circumference of the arm (measured in the center). The permissible range is printed on the cuff. If this is not suitable for your use, please call 1-866-464-6184. Large cuff sizes are sold separately.

Note: Only use clinically approved Microlife and CVS Cuffs!

- A loose cuff or a sideways protruding air pocket causes false measurement values.
- With repeated measurements, blood accumulates in the arm, which can lead to false results. Consecutive blood pressure measurements should be repeated after a 1 minute pause or after your arm has been held up in order to allow the accumulated blood to flow away.





Please note, clinical studies have proven the Microlife Averaging Mode (MAM) system with a 15 second wait period between measurements is a more reliable way to take your blood pressure. If you decide to retake your MAM measurement again, be sure to wait at least one minute before repeating a MAM measurement.

5.3. Fitting the cuff

- Pass the end of the cuff through the flat metal ring so that a loop is formed. The velcro closer must be facing outward. (Ignore this step if the cuff has already been prepared.)
- Place the cuff over the left upper arm so that the tube is closer to your lower arm.
- Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately $3/4$ " to 1" (2 to 3cm) above the elbow and that the tube is closer to the inner side of the arm.

Important! The **small white arrow on the cuff** must lie exactly over the artery which runs down the inner side of the arm.


- Tighten the cuff by pulling the end and close the cuff by affixing the velcro.
- There should be little free space between the arm and the cuff. You should be able to fit 2 fingers between your arm and the cuff. Clothing must not restrict the arm. Any piece of clothing which does must be removed. Cuffs that don't fit properly result in false measurement values. Measure your arm circumference if you are not sure of proper fit. Larger cuffs are available (Section 3B).
- Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked.
- Remain seated quietly for two minutes before you begin the measurement.

Comment:

If it is not possible to fit the cuff to your left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.

Comparable blood pressure measurements always require the same conditions (Relax for several minutes before a reading).


Cuff on right arm

5.4. Setting the time and date


This blood pressure monitor incorporates an integrated clock with date display. Time and date are recorded along with pressure and pulse. After new batteries have been inserted, the clock begins to run from the following setting: 2000-01-01 00:00 O'clock.

You must then re-enter the date and current time. Please proceed as follows (Example: Entering 2002-06-15 Time 09:30 AM):


1. Press the Time button for at least 3 seconds. Release the button and then press again. The display now indicates the set year, during which the four characters blink. If you insert new batteries the year may already be flashing.

2. The correct year can be entered by pressing the MEMORY button. Once for each change. (Example: 1 x press).


3. Press the Time button again. The display now switches to the current date, during which the first character (month) blinks.


Note: Holding the button down speeds up the procedure.


4. The current month can now be entered by pressing the MEMORY button. (Example: pressing 6 x advances 6 months)


5. Press the Time button again. The last two characters (day) are now blinking.

6. The current day can now be entered by pressing the MEMORY button. (Example: 15 x presses advances the day from the 1st to 15th)

7. Press the Time button again. The display now switches to the current time, during which the first character (hour) blinks.

9. Press the User/Time button again. The last two characters (Minutes) now blink.
10. The minutes can now be entered by pressing the MEMORY button. (Example: 30 x presses advances the time from 0 to 30 minutes)
11. Finally, press the time button once to set.

5.5. Reading the set date



After all settings have been made, click the User/Time button once. The date is briefly displayed and then the time. The input is now confirmed and the clock begins to run.

5.6. Selection of Microlife Average Mode

Clinical studies demonstrate taking multiple blood pressure readings and calculating an "average" is more likely to determine your true blood pressure. Your Microlife Averaging Mode unit allows you to switch the unit to a special Averaging Mode setting (MAM) that automatically takes multiple readings!

Microlife Averaging Mode / MAM:

- a) If you would like to take an Average Mode measurement, please slide the switch to the right toward MAM.
- b) Average Mode takes generally 3 measurements in succession and calculates the result.
- c) A  symbol in the display indicates that the unit is set to the Averaging Mode.

Single mode:

- a) If you would like to take single measurement, please slide the switch toward no. 1.
- b) Single mode only has 1 measurement

5.7. Selecting the User (User 1 or User 2)

1. Press and hold the TIME/DATE button for 3 seconds. The User 1 icon will flash in the lower left corner of the LCD screen.
2. Press the MEMORY button to toggle between User 1 and User 2.
3. Press TIME/DATE again to make your selection.

5.8. Measuring procedure

After the cuff has been appropriately positioned the measurement can begin:

- Press the ON/OFF/START button. The pump begins to inflate the cuff. In the display, the increasing cuff pressure is continually displayed.
- After automatically reaching an individual pressure, the pump stops and the pressure slowly falls. The cuff pressure is displayed during the measurement.
- When the device has detected your pulse, the heart symbol in the display begins to blink and a beep tone is audible for every pulse beat.
- When the measurement has been concluded, a long beep tone sounds. The measured systolic and diastolic blood pressure values, as well as the pulse are now displayed.
- The appearance of this symbol  signifies that an irregular heartbeat was detected. This indicator is only a caution. It is important that you be relaxed, remain still and do not talk during measurements. NOTE: We recommend contacting your physician if you see this  indicator frequently.
- The measurement results are displayed until you switch the device off. If no button is pressed for 5 minutes, the device switches off automatically.
- When the unit is set to the Microlife Average Mode setting, 3 separate measurements will take place in succession after which your result is calculated. There is a 15 second resting time in-between each measurement. A count down indicates the remaining time and a beep will sound 5 seconds before the 2nd and 3rd readings will begin.

Pumping
Pressure


Measuring

Systolic

Diastolic

Pulse


Measurement
complete

If one of the measurements causes an error message, it will be repeated one more time. If any additional error occurs, the measurement will be discontinued and error code displayed.

Expanding the MAM measurement to see the three individual readings - This function allows you to view the three individual measurements used to calculate the MAM measurement. After taking a MAM measurement hold the memory button for 3 seconds until you hear a short beep. (Do not hold the button longer than 7 seconds or you will delete all the readings in the memory.) Let go of the button and watch the screen. It will automatically scroll through the three measurements used in the MAM measurement.

5.9. Memory – displaying the last 99 measurements

The measured results are stored in the instrument until a new measurement is carried out. With the unit in the OFF position, press and hold the MEMORY button for at least 3-seconds. The screen will first show the average of all readings in memory (icon shown is "A") followed by the most recent readings along with time and date. The MR symbol is displayed in the lower left. Release the button. The last reading will be displayed. Each time you press the memory button an earlier measurement will be displayed along with time and date. If an irregular heartbeat interval was detected in any measurements stored in memory, the  icon will be displayed.

Before you delete all readings stored in the memory, make sure you will not need refer to the readings at a later date. Keeping a written record is prudent and may provide additional information for your doctor's visit.

In order to delete all stored readings, depress the MEMORY button for at least 7 seconds, the display will show the symbol "CL" and 3 short beep sounds will be heard to indicate deletion of stored readings.

5.10. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g the patient feels unwell), the ON/OFF button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.

6. Error Messages/Troubleshooting

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed. We recommend you switch the device off and take a new measurement making sure you allow time for relaxation. (example: error no. 1)

Error No.	Possible cause(s) / Solutions
ERR 1	The tube may have loosened, or no pulse was detected. *Ensure cuff connections are tight with proper cuff placement. See section (5.3)
ERR 2	Unnatural pressure impulses. Reason: The arm was moved during the measurement (Artefact). Repeat measurement keeping arm still.
ERR 3	Repeat measurement keeping arm still. If inflation of the cuff takes too long, the cuff is not correctly seated or the hose connection's not tight. Check connections and repeat.
ERR 5	The difference between systolic and diastolic is excessive. Measure again carefully following proper cuff procedures and ensure measurement under quiet conditions.
ERR 6	Single data differ too much during Average Mode even after 4 cycles. No average result can be displayed.

Other possible errors and their solutions

If problems occur when using the device, the following points should be checked:

Malfunction

Remedy

The display remains blank when the instrument is switched on although the batteries are in place.

- Check battery installation.
- If the display is unusual, remove the batteries and then exchange them for new ones. Check polarity.

The pressure does not rise although the pump is running.

- Check the connection of the cuff tube and connect properly.

The device frequently fails to measure, or the values measured are too low or high.

1. Fit the cuff correctly on the arm. (see 5.3)
2. Before starting measurement make sure that the cuff is not too tight and that clothing is not exerting pressure on the arm. Take articles of clothing off if necessary.
3. Measure blood pressure again in complete peace and quiet.

Every measurement results in different values, although the device functions normally and normal values are displayed.

- Please read the following information and the points listed under "Common sources of error" (5.2) Repeat the measurement.
- Blood pressure changes constantly. The observed readings may accurately reflect your pressure.

Blood pressure values differ from those measured by my doctor

- Record daily values and consult your doctor.
- Pressure readings in your doctor's office may be higher due to anxiety.

After the instrument has inflated the cuff the pressure falls very slowly, or not at all. (No reasonable measurement possible).

- Check cuff connections. Ensure the unit has not been tampered with.

8. Care and maintenance

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle this cuff carefully and avoid all types of stress through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. **The cuff must not be washed in a dishwasher, clothes washer, or submerged in water.**
- d) Handle the tube carefully. Do not pull on it. Do not allow the tubing to kink and keep it away from sharp edges.
- e) Do not drop the monitor or treat it roughly in any way. Avoid strong vibrations.
- f) **Never open the monitor!** This invalidates the manufacturer's warranty.

8. Warranty

Your CVS/Microlife blood pressure monitor is guaranteed **for a lifetime** against manufacturers' defects for the original purchaser only, from date of purchase. The warranty does not apply to damage caused by improper handling, accidents, professional use, not following the operating instructions or alterations made to the instrument by third parties.

There are no user serviceable parts inside. Batteries or damage from old batteries is not covered by the warranty. The warranty is only valid upon presentation of the warranty card. The cuff is warranted for a period of 1 year.

Please note: According to international standards, your monitor should be checked for calibration every 2 years

If fluctuations in readings are larger than 15 mmHg, and/or you hear irregular pulse tones, and/or the Irregular Heartbeat Detector icon appears frequently, consult your doctor.

In order to receive market clearance from governmental bodies, this device was subjected to strict clinical tests. The computer program used to measure blood pressure values was tested by experienced cardiac specialists in Germany.

The manufacture of your blood pressure monitor is in accordance with the terms of the European standard for blood pressure measuring devices (see technical data) under the supervision of the Technical Monitoring Association Essen (RWTüV-Essen).

Never attempt to repair the instrument yourself!

Any unauthorized opening of the instrument invalidates all warranty claims!

9. Certifications

Device standard:

This device is manufactured to meet the European and United States standards for non-invasive blood pressure monitors:

EN1060-1 / 1995

EN1060-3 / 1997

EN1060-4 / 2004

AAMI/ANSI SP10

Electromagnetic compatibility:

Device fulfills the stipulations of the International standard IEC60601-1-2

Clinical testing:

Clinical performance tests were carried out in the U.K. "Association Protocol and British HTP Protocol."

The B.H.S. (British Hypertension Society) has graded this product "Recommended for clinical and home use." This is the highest grading available for blood pressure monitors. Please see the B.H.S. web site. (Uses the same algorithm as B.H.S. graded "AA" model number BP3BT0-A)

10. Technical specifications:

Weight:	483 g (with batteries)
Size:	124 (W) x 205 (L) x 81 (H) mm
Storage temperature:	-20 to +50°C (23°F - 122°F)
Humidity:	15 to 90% relative humidity maximum
Operation temperature:	10 to 40°C (50°F - 104°F)
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	oscillometric
Pressure sensor:	capacitive
Measuring range:	
SYS/DIA:	30 to 280 mmHg
Pulse:	40 to 200 per minute
Cuff pressure display range:	0–299 mmHg
Memory:	Automatically stores the last 99 measurements.
Measuring resolution:	1 mmHg
Accuracy:	pressure within ± 3 mmHg pulse ± 5 % of the reading
Power source:	a) 4 AA batteries, 1.5 V b) AC adapter 6 V DC 600 mA (voltage 4.5 V DC to 6 V DC)
Accessories:	cuff type S101 for arm circumference 24 - 33.7 cm (9.5" - 13.25") (included) cuff type S102 for arm circumference 30.5–40.64 cm (12" - 16") (call Microlife for details)

Made in China

Technical alterations reserved!

US patents:
6,485,429 (Irregular Heartbeat Feature)
6,447,457 (Microlife Average Mode)

12. How to Contact Us

Distributed by:

CVS Pharmacy, Inc.

Woonsocket, RI 02895

©2004 CVS/pharmacy

Toll Free: 1-866-464-6184

Email Questions to: CVSbpsupport@microlifeusa.com

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Manufacturado por Microlife Corporation