AUTOMATIC INFLATE Blood Pressure Monitor

Product #271245
Instruction Manual
Instrucciones de Uso

Questions? Preguntas?
1-866-464-6184

PLEASE NOTE:
THIS MEDICAL INSTRUMENT MUST BE USED ACCORDING
TO INSTRUCTIONS TO ENSURE ACCURATE READINGS.

CVS Automatic Blood Pressure Monitor

Instruction Manual

Table of contents

1. Introduction

- 1.1. Features of your blood pressure monitor
- 1.2. Important information about self-measurement

2. Important information on blood pressure and its measurement

- 2.1. How does high/low blood pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done if regular high/low values are obtained?

3. Components of your blood pressure monitor

4. Operation of your blood pressure monitor

- 4.1. Inserting the batteries
- 4.2. Using an A/C power adapter
- 4.3. Tube connection
- 4.4. Setting the Time and Date

5. Measurement Procedure

- 5.1. Before measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measuring procedure
- 5.5. Reading Recall – displaying the last measurement
- 5.6. Discontinuing a measurement

6. Error Messages/Troubleshooting

7. Care and Maintenance

8. Warranty

9. Certifications

10. Technical specifications

11. How to contact us

1. Introduction

1.1. Features of your CVS automatic blood pressure monitor, model #271245

Your blood pressure monitor is a fully automatic, digital blood pressure measuring device for use on the upper arm, which enables very fast and reliable measurement of the systolic and diastolic blood pressure as well as the pulse by way of the oscillometric method. This device offers clinically proven accuracy and has been designed to be user friendly.

Before using, please read this instruction manual carefully and then keep it in a safe place. Please contact your doctor for further questions on the subject of blood pressure measurement.

Attention!

1.2. Important information about self-measurement

- **Self-measurement means Control**, not diagnosis or treatment. Unusual values must always be discussed with your doctor. **Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.**
- The pulse display is **not** suitable for checking the frequency of heart pacemakers!
- In cases of irregular heartbeat (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with your doctor.

Electromagnetic interference:

The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring accuracy.

2. Important information on blood pressure and its measurement

2.1. How does high/low blood pressure arise?

Your level of blood pressure is determined in the circulatory center of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is effected by fine muscles in the blood vessel walls.

Your level of arterial blood pressure changes periodically during heart activity: During the “blood ejection” (Systole) the value is highest (systolic blood pressure value). At the end of the heart’s “rest period” (Diastole) pressure is lowest (diastolic blood pressure value). Blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2. Which values are normal?

Blood pressure is too high **if at rest**, your diastolic pressure is above 90mmHg and/or the systolic blood pressure is over 160mmHg. In this case, please consult your doctor immediately. Long-term values at this level endanger your health due to continual damage to the blood vessels in your body.

Should the systolic blood pressure values lie between 140mmHg and 160mmHg and/or the diastolic blood pressure values lie between 90mmHg and 95mmHg. Consult your doctor. Regular self-checks will be necessary.

With blood pressure values that are too low, (i.e. systolic values under 105mmHg and/or diastolic values under 60mmHg) consult your doctor.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your doctor. **Never use the results of your measurements to independently alter the drug doses prescribed by your doctor.**

The following standards for assessing high blood pressure (without regard to age) have been established by the National Institutes of Health JNCVI.

Category	Systolic (mmHg)	Diastolic (mmHg)
Optimal	<120	<80
Normal	<130	<85
High Normal	130-139	85-89
Hypertension		
Stage 1	140-159	90-99
Stage 2	160-179	100-109
Stage 3	≥180	≥110

☞ **Further information**

- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called “labile hypertension”. Consult your doctor.
- **Correctly measured diastolic blood pressure values above 120mmHg require immediate medical treatment.**

2.3. What can be done if regular high or low values are obtained?

- a) Consult your doctor.
- b) Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result from arteriosclerosis. Furthermore, the heart will become structurally damaged with increased diastolic blood pressure values.
- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organic malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take to reduce and prevent high blood pressure. These measures must be permanent lifestyle changes.

A) Eating habits

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

B) Previous illnesses

- Consistently follow all medical instructions for treating illnesses such as:
 - Diabetes (Diabetes mellitus or sugar diabetes)
 - Fat metabolism disorder
 - Gout

C) Habits

- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (coffee, tea, chocolate)


D) Physical constitution

- After a preliminary medical examination, exercise regularly.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your sporting activities. He/She will advise you regarding the type and extent of types of sport that are possible for you.

3. Components of your blood pressure monitor

The illustration shows the blood pressure monitor consisting of:


a) Measuring unit

b) Upper arm cuffs:


The cuff is for arm circumference 24–34 cm or
9.5" - 13.25" (included)

CVS sells a large cuff for arm circumference 31–41 cm or
12" - 16" (sold separately in all CVS stores)


Please Note: Do not force cuff connection into the opening. Make sure the cuff connection is not pushed into the AC adapter port. If the cuff is too small, call toll free 1-866-464-6184 for further information.

Please note: The AC adapter port is in the back of the unit and the cuff port is on the side. See illustrations to the right.

Cuff port


AC Adapter port

4. Operation of your blood pressure monitor

4.1. Inserting the batteries

After you have unpacked your device insert the batteries. The battery compartment is located on the back side of the device (see illustration).

- Remove cover as illustrated
- Insert the batteries (4 x size AA 1,5 V), following the indicated polarity. (+ -)
- If the low battery indicator appears in the display, the batteries are discharged and the unit will not function.
- If the batteries are inserted incorrectly, the display may function erratically or not at all. Check battery polarity.

Attention!

- After the low battery indicator appears, the device won't function until the batteries have been replaced.
- Please use "AA" Long-Life or Alkaline 1.5V Batteries.
- Do not use rechargeable batteries.
- If the blood pressure monitor is not used for long periods, remove the batteries from the device.

Functional check: Hold on the On/Off button down to test all the display elements. When functioning correctly all icons will appear.


Low Battery Indicator

4.2. Using an AC power adaptor (special accessory sold separately. Call 1-866-464-6184 to order)

It is possible to operate this blood pressure instrument with an AC adaptor (output 6 V DC / 600 mA with DIN plug). Make certain that you use an adaptor which fulfils the legal requirements and electronic requirements in the U.S. (UL standard)

- Push the plug into the socket on the back side of the instrument.
- Push the AC adaptor into a 110 V power socket (U.S. or Canada).
- Test that power is available by pressing the On/Off button.


Note:

- No power is taken from the batteries while the AC adaptor is connected to the instrument.
- If the power is interrupted during the measurement (e.g. by removal of the adaptor from the wall socket), the instrument must be reset by removing the plug from the instrument.
- Please consult CVS if you have questions relating to the AC adaptor.

4.3. Tube connection

Insert the cuff tube into the opening provided on the left side of the instrument, as shown in the diagram.

4.4. Setting the Time and Date

To set the time and date press the Set Time/Date button and hold it down for five seconds. When you let go the number you'll see a flashing four digit number (probably 2001 or 2002 to 2006). This number represents the year. Use the MEMORY key to scroll up (you cannot scroll backwards so scroll slowly) to the correct year. When you reach the correct year stop scrolling and press the Set Time/Date button to set the year.

After doing this the number flashing will be the month. Again, use the MEMORY key to scroll to the appropriate month and stop scrolling then press the Set Time/Date button to set the month and move onto setting the day. Use the MEMORY key to scroll up to find the correct day of the month and press the Set Time/Date button to set it.

Now that you've set the correct date, the number flashing will represent the time. The first two digits are the hour. Use the MEMORY button to scroll through the hours. Note that for morning hours before noon, there is no AM displayed but after noon PM will display. Again, scroll through the time using the MEMORY button to scroll and the Set Time/Date button to set the time.

5. Measurement Procedure

Please note: You should always be seated before and during measurement

5.1. Before the measurement:

- Avoid eating and smoking as well as all forms of exertion directly before measurement. These factors influence the measurement result. Find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Always perform measurements at the same time of day, since the blood pressure changes during the course of the day.

5.2. Common sources of error:

Note:



Comparable blood pressure measurements always require the same conditions!

These are normally quiet conditions.



- Efforts by the patient to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support and rest your arm on a table.
- If the upper arm artery lies considerably lower or higher than the heart, a false high or low blood pressure will be measured!
***A variation of 6" between the cuff and heart level can result in a reading error of + or - 10mmHg.**
- Cuffs that do not fit properly result in false measurement values. Selecting the correct cuff is extremely important. The cuff size is dependent upon the circumference or distance around your upper arm measured in the center. The permissible range is printed on the cuff. If this is not suitable for your use, please contact CVS.
Note: Only use clinically approved **Microlife or CVS Cuff Model # 271237.**
- A loose cuff or a sideways protruding air-pocket causes false measurement values.
- Measurements should be done after a 5 minute rest to ensure accuracy.
- Consecutive blood pressure measurements should be performed after a 1 minute pause.

5.3. Fitting the cuff

- Pass the end of the cuff through the flat metal ring so that a loop is formed. The Velcro closer must be facing outwards. (Ignore this step if the cuff has already been prepared.)
- Place the cuff over the left upper arm so that the tube is closer to your lower arm. Align the cuff artery mark (arrow) over the arm artery than runs down the inner side of the arm.
- Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately 3/4" to 1" (2 to 3 cm) above the elbow and that the tube is closer to the inner side of the arm.

- d) Tighten the cuff by pulling the end and close the cuff by affixing the velcro. Make sure that the cuff Index end falls between the OK lines when the cuff is in position. If it doesn't, you are using the wrong cuff size.
- e) There should be little free space between the arm and the cuff. You should be able to fit 2 fingers between your arm and the cuff. Clothing must not restrict the arm. Any piece of clothing which does must be removed.
- f) Lay your arm on a table (palm upward) so the cuff is at the same height as the heart. Make sure that the tube is not kinked.
- g) Remain seated quietly for five minutes before you begin the measurement.

Cuff on right arm

Comment:

If it is not possible to fit the cuff to the left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.


5.4. Measuring procedure

After the cuff has been appropriately positioned the measurement can begin:

- a) Press the On/Off/Memory button. The pump begins to inflate the cuff. In the display, the increasing cuff-pressure is continually displayed.
- b) After automatically reaching an individual pressure, the pump stops and the pressure slowly falls. The cuff-pressure is displayed during the measurement.



- c) When the device has detected your pulse, the heart symbol in the display begins to blink and a beep tone is audible for every pulse beat.
- d) When the measurement has been concluded, a long beep tone sounds. The measured systolic and diastolic blood pressure values as well as the pulse are now displayed.
- e) The measurement results are displayed until you switch the device off. If no button is pressed for 5 minutes, the device switches off automatically.


5.5. Memory – displaying the last measurement

The measured results are stored in the instrument. This unit stores 99 measurements. You can call up the stored values at any time by pressing the MEMORY button. The first memory displayed will show an A for average. This is an average of all measurements stored in the memory. If you press the memory a second time, you'll see the most recent measurement taken. Keep pushing the button to scroll back through memories.

5.6. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g the patient feels unwell), the On/Off button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.


Further information

Measurements should not occur soon after each other, since the results will be inaccurate. Wait several minutes in a relaxed position, sitting or lying, before you repeat a measurement.

6. Error Messages/Troubleshooting

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed.
(Example Error no. 1)

Error No.	Possible cause(s)
ER 1	The systolic pressure was determined but afterwards the cuff pressure fell below 20 mmHg. This situation can occur, for example, if the tube becomes unplugged after the systolic blood pressure has been measured. Further possible cause: No pulse has been detected.
ER 2	Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the measurement (Artefact).
ER 3	Inflation of the cuff takes too long. The cuff is not correctly seated or the hose connection is not sealed.
ER 5	The difference between systolic and diastolic is excessive. Measure again carefully following proper cuff procedures and ensure measurement under quiet conditions. Contact your doctor if this message occurs frequently.

Other possible errors and their solutions

If problems occur when using the device, the following points should be checked:

Malfunction	Remedy
The display remains blank when the instrument is switched ON although the batteries are in place.	<ul style="list-style-type: none">• Check battery installation.• If the display is unusual, remove the batteries and then exchange them for new ones• Check polarity
The pressure does not rise although the pump is running.	<ul style="list-style-type: none">• Check the connection of the cuff tube and connect properly.
The device frequently fails to measure the blood pressure values, or the values measured are too low or too high.	<ol style="list-style-type: none">1. Fit the cuff correctly on the arm. (see 5.3)2. Before starting the measurement make sure that the cuff is not too tight and that clothing is not exerting pressure on the arm. Take off articles of clothing if necessary.3. Measure blood pressure again in complete peace and quiet.
Every measurement results in different values, although the device functions normally and normal values displayed.	<ul style="list-style-type: none">• Please read the following information and the points listed under "Common sources of error"(5.2). Repeat the measurement.• Blood Pressure changes constantly. The observed readings may accurately reflect your pressure.
The blood pressure values differ from those measured by my doctor.	<ul style="list-style-type: none">• Record daily values and consult your doctor.• Pressure readings in your doctor's office may be higher due to anxiety.• Your blood pressure fluctuates throughout the day so record the time along with your reading.
After the instrument has inflated the cuff the pressure falls very slowly or not at all. (No reasonable measurement possible).	<ul style="list-style-type: none">• Check cuff connections. Ensure the unit has not been tampered with.

7. Care and maintenance

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds.
The cuff must not be washed in a dishwasher, clothes washer, or submerged in water.
- d) Handle the tube carefully. Do not pull on it. Do not allow the tubing to kink and keep it away from sharp edges.
- e) Do not drop the monitor or treat it roughly in any way. Avoid strong vibrations.
- f) **Never open the monitor!** This invalidates the manufacturer's warranty.

Further information

Blood pressure is subject to fluctuations even in healthy people.

Comparable measurements always require the same conditions (quiet conditions)!

If fluctuations in readings are larger than 15mmHg, and/or you hear irregular pulse tones, consult your doctor.

In order to receive market clearance from governmental bodies, this device was subjected to strict clinical tests. The computer program used to measure blood pressure values was tested by experienced cardiac specialists in Germany.

The manufacture of your CVS blood pressure monitor is in accordance with the terms of the European standard for blood pressure measuring devices (see technical data) under the supervision of the Technical Monitoring Association Esses (RWTüV-Essen).

Never attempt to repair the instrument yourself!

Any unauthorized opening of the instrument invalidates all warranty claims!

8. Warranty

Your CVS blood pressure monitor is **guaranteed for a lifetime** against manufacturers' defects for the original purchaser only, from date of purchase. The warranty does not apply to damage caused by improper handling, accidents, professional use, not following the operating instructions or alterations made to the instrument by third parties.

There are no user serviceable parts inside. Batteries or damage from old batteries is not covered by the warranty. The warranty is only valid upon completion and delivery of the warranty card to address on the enclosed warranty card.

Please call toll free 1-866-464-6184 or email CVSbpsupport@microlifeusa.com for warranty claims.

Please note: According to international standards, your monitor should be checked for calibration every 2 years.

9. Certifications

Device standard:

This device is manufactured to meet the European and United States standards for non-invasive blood pressure monitors
EN1060-1 / 1995
EN1060-3 / 1997

DIN 58130, NIBP – clinical investigation
ANSI / AAMI SP10, NIBP – requirements

Electromagnetic compatibility:

Device fulfills the stipulations of the European standard EN 60601-1-2

Clinical testing:

Clinical performance tests were carried out in the US and Germany according to the DIN 58130/1997 procedure N6 (sequential) and AAMI standard (US).

The stipulations of the EU-Guidelines 93/42/EWG for Medical Products Class IIa have been fulfilled.

10. Technical specifications:

Weight:	465 g (with batteries)
Size:	131 (W) x 174 (L) x 73 (H) mm
Storage temperature:	-5 to +50°C
Humidity:	15 to 85% relative humidity maximum
Operation temperature:	10 to 40°C
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	oscillometric
Pressure sensor:	capacitive
Measuring range:	
SYS/DIA:	30 to 280 mmHg
Pulse:	40 to 200 per minute
Cuff pressure display range:	0-299 mmHg
Memory:	Storing the last 99 measurements automatically
Measuring resolution:	1 mmHg
Accuracy:	Pressure within ± 3 mmHg Pulse ± 5 % of the reading
Power source:	a) 4 dry cells (batteries) UM-3, size AA, 1.5 V b) AC adaptor 6 V DC 600 mA (voltage 4.5 V DC to 6 V DC)
Accessories:	regular cuff for arm circumference 9.5"-13.25" large cuff for arm circumference 12"-16"

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©2004 CVS/pharmacy
CVS Pharmacy, Inc.