microlife Digital Fully Automatic Blood Pressure Monitor with Patented Irregular Heartbeat Detector

Model # 3BT0-AP Instruction Manual

CAUTION: If the Heartbeat Interval Irregularity Indicator icon comes on frequently, see your doctor. You may require medical attention. Only a physician can diagnose irregular heartbeat.

US Patent # 6,485,429
(Irregular Heartbeat Detector)

microlife

Automatic Blood Pressure Monitor

Instruction Manual

Table of contents

1. Introduction

- 1.1. Features of your blood pressure monitor
- 1.2. Important information about self measurement

2. Important information on blood pressure and its measurement

- 2.1. How does high/low blood pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done if regular high/low values are obtained?
- 2.4. Irregular Heartbeat Detector

3. Components of your blood pressure monitor

4. Operation of your blood pressure monitor

- 4.1. Inserting the batteries
- 4.2. Using an A/C power adapter
- 4.3. Tube connection

Measurement Procedure

- 5.1. Before measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Setting the time and date
- 5.5. Measuring procedure
- 5.6. Memory – displaying the last measurement
- 5.7. Discontinuing a measurement

6. Error Messages/Troubleshooting

7. Care and maintenance

8. Warranty

9. Certifications

10. Technical specifications

11. How to Contact Us

1. Introduction

1.1. Features of your Microlife automatic blood pressure monitor, model # 3BT0-AP

Your Microlife blood pressure monitor is a fully automatic digital blood pressure measuring device with an integrated Irregular Heartbeat Detector, for use by adults on the upper arm at home or in your doctor's/nurse's office. It enables very fast and reliable measurement of the systolic and diastolic blood pressure as well as the pulse by way of the oscillometric method. This device offers clinically proven accuracy and has been designed to be user friendly.

Before using, please read this instruction manual carefully and then keep it in a safe place. Please contact your doctor for further questions on the subject of blood pressure and its measurement.

Attention!

1.2. Important information about self-measurement

- **Self-measurement means Control**, not diagnosis or treatment. Unusual values must always be discussed with your doctor. **Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.**
- The pulse display is **not** suitable for checking the frequency of heart pacemakers!
- In cases of irregular heartbeat (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with your doctor.

Electromagnetic interference:

The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring accuracy.

2. Important information on blood pressure and its measurement

2.1. How does high/low blood pressure arise?

Your level of blood pressure is determined in the circulatory center of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is effected by fine muscles in the blood vessel walls.

Your level of arterial blood pressure changes periodically during heart activity: During the "blood ejection" (Systole) the value is highest (systolic blood pressure value). At the end

of the heart's "rest period" (Diastole) pressure is lowest (diastolic blood pressure value). Blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2. Which values are normal?

Blood pressure is too high if your diastolic pressure is above 90 mmHg and/or your systolic blood pressure is over 160 mmHg, **while at rest**. In this case, please consult your physician immediately. Long-term values at this level endanger your health due to continual damage to the blood vessels in your body.

If your systolic blood pressure values are between 140 mmHg and 160mmHg and/or the diastolic blood pressure values between 90 mmHg and 95mmHg, consult your physician. Regular self-checks are necessary.

If you have blood pressure values that are too low, (i.e., systolic values under 105mmHg and/or diastolic values under 60 mmHg), consult your physician.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your physician. **Never use the results of your measurements to independently alter the drug doses prescribed by your physician.**

Which values are normal?

The following standards for assessing high blood pressure (without regard to age) have been established by the National Institutes of Health JNCVI.

Category	Systolic (mmHg)	Diastolic (mmHg)
Optimal	<120	<80
Normal	<130	<85
High Normal	130-139	85-89
Hypertension		
Stage 1	140-159	90-99
Stage 2	160-179	100-109
Stage 3	≥180	≥110

 **Further information**

- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called “labile hypertension.” Consult your doctor.
- **Correctly measured diastolic blood pressure values above 120mmHg and require immediate medical treatment.**

2.3. What can be done if regular high or low values are obtained?

- a) Consult your doctor.
- b) Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result from arteriosclerosis. Furthermore, the heart will become structurally damaged with increased blood pressure values.
- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organ malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take to reduce and even prevent high blood pressure. These measures must be permanent lifestyle changes.

A) Eating habits

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

B) Previous illnesses

- Consistently follow all medical instructions for treating illness such as:
 - Diabetes (Diabetes mellitus or sugar diabetes)
 - Fat metabolism disorder
 - Gout

C) Habits

- Give up smoking completely.
- Drink only moderate amounts of alcohol.
- Restrict your caffeine consumption (coffee, tea, chocolate).

D) Your Physical condition

- **After a medical examination, and with your doctor's approval and direction; exercise.**
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your exercise routine. You must receive advise regarding the type and extent of exercise that is appropriate for you.

2.4. Irregular Heartbeat Detector


The appearance of this symbol  signifies that a certain irregular heartbeat was detected which could cause cardiovascular disease. This indicator is only a caution. The Microlife Fully Automatic Blood Pressure Monitor with Heartbeat Irregular Heartbeat Detector can detect blood pressure and pulse rate measurement with specifications even when the patient is having an irregular heartbeat interval. It is important that you be relaxed, remain still and do not talk during the measurements.

Note: We recommend contacting your physician if you see this '' icon frequently.

CAUTION: An irregular heartbeat can be a serious medical condition requiring medical attention. See your doctor if you suspect you have an irregular heartbeat.

3. Components of your blood pressure monitor

a) Measuring unit

b) Upper arm cuffs:

Type S101 for arm circumference 22-33.7 cm or

9.5" - 13.25" (included)

Type S102 for arm circumference 30.5-40.64 cm

12" - 16" (order as a special accessory)


Please Note:

Arm circumference should be measured with a measuring tape in the middle of the relaxed upper arm. Do not force cuff connection into the opening. Make sure the cuff connection is not pushed into the AC adapter port. If the cuff is too small, call *microlife* for further information.

4. Operation of your blood pressure monitor**4.1. Inserting the batteries**

After you have unpacked your device insert the batteries. The battery compartment is located on the back side of the device (see illustration).

- a) Remove cover as illustrated.
- b) Insert the batteries (4 AA, 1.5 V), following the indicated polarity. (+ -)
- c) If the low battery indicator appears in the display, the batteries are discharged and the unit will not function.
- d) If batteries are inserted incorrectly, the display may function erratically or not at all. Check battery polarity.

Low Battery Indicator

Attention!


- After the low battery indicator appears, the device won't function until the batteries have been replaced.
- Please use "AA" Long-Life or Alkaline 1.5V Batteries.
- Do not use rechargeable batteries.
- If the blood pressure monitor is not used for long periods, remove the batteries from the device.

Functional check: Hold the On/Off button down to test all the display elements. When functioning correctly many icons will appear.

4.2. Using an AC power adapter (special accessory)

It is possible to operate this blood pressure instrument with an AC adapter. (output 6 V DC / 600 mA with DIN plug). Make certain that you use an adapter which fulfills the legal requirements and electronic requirements in the U.S. (UL standard). Use of an AC adapter other than a Microlife adapter may invalidate your warranty.


- Push the plug into the socket at the left side of the instrument.
- Plug the AC adapter into a 110 V power socket (U.S. or Canada).
- Test that power is available by pressing the On/Off button.

**Note:**

- No power is taken from the batteries while the AC adapter is connected to the instrument.
- If the power is interrupted during a measurement (e.g. by removal of the adapter from the wall socket), the instrument must be reset by removing the plug from the instrument.
- Please call 1-800-568-4147 if you have questions relating to the AC adapter.

4.3. Tube connection

Insert the cuff tube into the opening provided on the side of the instrument, as shown in the diagram.

5. Measurement Procedure

Please note: You should always be seated before and during measurement.

5.1. Before measurement:

- Avoid eating and smoking as well as all forms of exertion directly before measurement. These factors influence the measurement result. Find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Always perform measurements at the same time of day, since blood pressure changes during the course of the day.

5.2. Common sources of error:

Note: Comparable blood pressure measurements always require the same conditions! Conditions should always be quiet.



- All efforts by the user to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not flex any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower or higher than the heart, an erroneously high or low blood pressure will be measured! Each 15cm (6 inch) difference in height between your heart and the cuff results in a measurement error of 10 mmHg!
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is extremely important. Cuff size is dependent upon the circumference of the arm (measured in the center). The permissible range is printed on the cuff. If this is not suitable for your use, please contact microlife USA. Large cuff sizes are sold separately.

Note: Only use clinically approved microlife Cuffs!

- A loose cuff or a sideways protruding air pocket causes false measurement values.
- With repeated measurements, blood accumulates in the arm, which can lead to false results. Consecutive blood pressure measurements should be repeated after a 1 minute pause or after your arm has been held up in order to allow the accumulated blood to flow away.




5.3. Fitting the cuff

- Pass the end of the cuff through the flat metal ring so that a loop is formed. The velcro closer must be facing outward. (Ignore this step if the cuff has already been prepared.)
- Place the cuff over the left upper arm so that the tube is closer to your lower arm.
- Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately $\frac{3}{4}$ " to 1" (2 to 3cm) above the elbow and that the tube is closer to the inner side of the arm.

Important! The **small white arrow on the cuff** must lie exactly over the artery which runs down the inner side of the arm.


- Tighten the cuff by pulling the end and close the cuff by affixing the velcro.
- There should be little free space between the arm and the cuff. You should be able to fit 2 fingers between your arm and the cuff. Clothing must not restrict the arm. Any piece of clothing which does must be removed. Cuffs that don't fit properly result in false measurement values. Measure your arm circumference if you are not sure of proper fit. Larger cuffs are available (Section 3B).
- Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked.
- Remain seated quietly for two minutes before you begin the measurement.

Comment:

If it is not possible to fit the cuff to your left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.

Comparable blood pressure measurements always require the same conditions (Relax for several minutes before a reading).

Cuff on right arm


5.4. Setting the time and date

This blood pressure monitor incorporates an integrated clock with date display. Time and date are recorded along with pressure and pulse. After new batteries have been inserted, the clock begins to run from the following setting: 2000-01-01 00:00 O'clock.


You must then re-enter the date and current time. Please proceed as follows (Example:

Entering 2002-06-15 Time 09:30 AM):


1. Press the Time button for at least 3 seconds. The display now indicates the set year, during which the four characters blink.

2. The correct year can be entered by pressing the MEMORY button. Once for each change. (Example: 1 x press).

3. Press the Time button again. The display now switches to the current date, during which the first character (month) blinks.


Note: Holding the button down speeds up the procedure.


4. The current month can now be entered by pressing the MEMORY button. (Example: pressing 6 x advances 6 months)


5. Press the Time button again. The last two characters (day) are now blinking.

6. The current day can now be entered by pressing the MEMORY button. (Example: 15 x presses advances the day from the 1st to 15th)

7. Press the Time button again. The display now switches to the current time, during which the first character (hour) blinks.

9. Press the User/Time button again. The last two characters (Minutes) now blink.
10. The minutes can now be entered by pressing the MEMORY button. (Example: 30 x presses advances the time from 0 to 30 minutes)

5.5. Reading the set date


After all settings have been made, click the User/Time button once. The date is briefly displayed and then the time. The input is now confirmed and the clock begins to run.

5.6. Measuring procedure


After the cuff has been appropriately positioned the measurement can begin:

- Press the ON/OFF/START button. The pump begins to inflate the cuff. In the display, the increasing cuff pressure is continually displayed.
- After automatically reaching an individual pressure, the pump stops and the pressure slowly falls. The cuff pressure is displayed during the measurement.
- When the device has detected your pulse, the heart symbol in the display begins to blink and a beep tone is audible for every pulse beat.
- When the measurement has been concluded, a long beep tone sounds. The measured systolic and diastolic blood pressure values, as well as the pulse are now displayed.
- The appearance of this symbol  signifies that a certain irregular heartbeat was detected. This indicator is only a caution. It is important that you be relaxed, remain still and do not talk during measurements. NOTE: We recommend contacting your physician if you see this  icon frequently.
- The measurement results are displayed until you switch the device off. If no button is pressed for 5 minutes, the device switches off automatically.


Pumping
Pressure

Measuring

Systolic
Diastolic
Pulse


Measurement
complete

5.7. Memory – displaying the last 99 measurements

The measured results are stored in the instrument until a new measurement is carried out. With the unit in the OFF position, press and hold the MEMORY button for at least 3-seconds. The screen will show the last reading along with time and date.

The **MR** symbol is displayed in the lower left. Release the button. The last reading will be displayed. Each time you press the memory button an earlier measurement will be displayed along with time and date. If an irregular heartbeat was detected in any measurements stored in memory, the  icon will be displayed.


5.8. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g the patient feels unwell), the ON/OFF button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.

6. Error Messages/Troubleshooting

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed. We recommend you switch the device off and take a new measurement making sure you allow time for relaxation.
(example: error no. 1)

Error No.	Possible cause(s) / Solutions
ERR 1	The tube may have loosened, or no pulse was detected. *Ensure cuff connections are tight with proper cuff placement. See section (5.3)
ERR 2	Unnatural pressure impulses. Reason: The arm was moved during the measurement (Artefact). Repeat measurement keeping arm still.

ERR 3	Repeat measurement keeping arm still. If inflation of the cuff takes too long, the cuff is not correctly seated or the hose connection's not tight. Check connections and repeat.
ERR 5	The difference between systolic and diastolic is excessive. Measure again carefully following proper cuff procedures and ensure measurement under quiet conditions.

Other possible errors and their solutions

If problems occur when using the device, the following points should be checked:

Malfunction	Remedy
The display remains blank when the instrument is switched on although the batteries are in place.	<ul style="list-style-type: none"> • Check battery installation. • If the display is unusual, remove the batteries and then exchange them for new ones. Check polarity.
The pressure does not rise although the pump is running.	<ul style="list-style-type: none"> • Check the connection of the cuff tube and connect properly.
The device frequently fails to measure, or the values measured are too low or high.	<ol style="list-style-type: none"> 1. Fit the cuff correctly on the arm. (see 5.3) 2. Before starting measurement make sure that the cuff is not too tight and that clothing is not exerting pressure on the arm. Take articles of clothing off if necessary. 3. Measure blood pressure again in complete peace and quiet.
Every measurement results in different values, although the device functions normally and normal values are displayed.	<ul style="list-style-type: none"> • Please read the following information and the points listed under "Common sources of error".(5.2) Repeat the measurement. • Blood pressure changes constantly. The observed readings may accurately reflect your pressure.
Blood pressure values differ from those measured by my doctor	<ul style="list-style-type: none"> • Record daily values and consult your doctor. • Pressure readings in your doctor's office may be higher due to anxiety.
After the instrument has inflated the cuff the pressure falls very slowly, or not at all. (No reasonable measurement possible).	<ul style="list-style-type: none"> • Check cuff connections. <p>Ensure the unit has not been tampered with.</p>

7. Care and maintenance

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle this cuff carefully and avoid all types of stress through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. **The cuff must not be washed in a dishwasher, clothes washer, or submerged in water.**
- d) Handle the tube carefully. Do not pull on it. Do not allow the tubing to kink and keep it away from sharp edges.
- e) Do not drop the monitor or treat it roughly in any way. Avoid strong vibrations.
- f) **Never open the monitor!** This invalidates the manufacturer's warrantee.

9. Warranty

Your microlife blood pressure monitor is **guaranteed for 5 years** against manufacturers' defects for the original purchaser only, from date of purchase. The warranty does not apply to damage caused by improper handling, accidents, professional use, not following the operating instructions or alterations made to the instrument by third parties.

There are no user serviceable parts inside. Batteries or damage from old batteries is not covered by the warranty. The warranty is only valid upon presentation of the warranty card.

Please note: According to international standards, your monitor should be checked for calibration every 2 years.

If fluctuations in readings are larger than 15 mmHg, and/or you hear irregular pulse tones, and/or the Heartbeat Interval Irregularity Indicator appears frequently, consult your doctor.

In order to receive market clearance from governmental bodies, this device was subjected to strict clinical tests. The computer program used to measure blood pressure values was tested by experienced cardiac specialists in Germany.

The manufacture of your *microlife* blood pressure monitor is in accordance with the terms of the European standard for blood pressure measuring devices (see technical data) under the supervision of the Technical Monitoring Association Essen (RWTüV-Essen).

Never attempt to repair the instrument yourself!

Any unauthorized opening of the instrument invalidates all warrantee claims!

9. Certifications

Device standard:

This device is manufactured to meet the European and United States standards for non-invasive blood pressure monitors:

EN1060-1 / 1995

EN1060-3 / 1997

ANSI / AAMI SP10, NIBP – requirements

Electromagnetic compatibility:

Device fulfills the stipulations of the International standard IEC60601-1-2

Clinical testing:

Clinical performance tests were successfully completed against AAMI SP10-1992 and British Hypertension Society protocol B.H.S.: 1993. Studies available on request.

The B.H.S. (British Hypertension Society) has rated this product "Recommended for clinical and home use." This is the highest grading available for blood pressure monitors. Please see the B.H.S. web site. (Uses the same algorithm as B.H.S. graded "AA" model number 3BT0-1)

10. Technical specifications:

Weight:	483 g (with batteries)
Size:	124 (W) x 205 (L) x 81 (H) mm
Storage temperature:	-5 to +50°C (23°F - 122°F)
Humidity:	15 to 85% relative humidity maximum
Operation temperature:	10 to 40°C (50°F - 104°F)
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	oscillometric
Pressure sensor:	capacitive
Measuring range:	
SYS/DIA:	30 to 280 mmHg
Pulse:	40 to 200 per minute
Cuff pressure display range:	0-299 mmHg
Memory:	Automatically stores the last 99 measurements.
Measuring resolution:	1 mmHg
Accuracy:	pressure within ± 3 mmHg pulse ± 5 % of the reading
Power source:	a) 4 AA batteries, 1.5 V b) AC adapter 6 V DC 600 mA (voltage 4.5 V DC to 6 V DC)
Accessories:	cuff type S101 for arm circumference 24-33.7 cm (9.5" -13.25") (included) cuff type S102 for arm circumference 30.5-40.64 cm (12" - 16") (call Microlife for details)

Made in China.

Technical alterations reserved!

Manufactured by Microlife Corporation

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Manufacturado por Microlife Corporation
Microlife USA, Inc.