**microlife Premier Wrist Blood Pressure Monitor**

Model# 3BU1-5

Instruction Manual

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Manual de instrucciones

PLEASE NOTE:  
THIS MEDICAL INSTRUMENT MUST BE USED ACCORDING TO  
INSTRUCTIONS TO ENSURE ACCURATE READINGS

Questions ? Call microlife toll-free at 1-800-568-4147

[www.microlifeusa.com](http://www.microlifeusa.com)

**CE0044**

## **Premier Wrist Blood Pressure Monitor**

Instruction Manual

Table of contents

---

### **1. Introduction**

---

- 1.1.Features of the BP 3BU1-5
- 1.2.Important information about self-measurement

### **2. Important information on the subject of blood pressure and its measurement**

---

- 2.1.How does high/low blood pressure arise?
- 2.2.Which values are normal?
- 2.3.What can be done, if regular high/low values are obtained?

### **3. The various components of the blood pressure monitor**

---

#### **4. Putting the blood pressure monitor into operation**

---

- 4.1.Inserting the batteries
- 4.2.Sleep mode
- 4.3.Setting the time and date
- 4.4.User selection

#### **5. Carrying out a measurement**

---

- 5.1.Before the measurement
- 5.2.Common sources of error
- 5.3.Fitting the cuff
- 5.4.Measuring procedure
- 5.5.Discontinuing a measurement
- 5.6.Memory – storage and recall of the measurements
- 5.7.Memory – recall of 3 intermittent average results
- 5.8.Memory – cancellation of all measurements

## **6. Printer functions**

- 6.1. How to print the last stored measurement (single data print)
- 6.2. How to print all measurements stored in memory (data print out with chart)
- 6.3. Stop printing

## **7. Error messages/malfunctions**

## **8. Care and maintenance**

## **9. Warranty**

## **10. Reference to standards**

## **11. Technical specifications**

## **12. [www.microlifeusa.com](http://www.microlifeusa.com)**

## **1. Introduction**

### **1.1. Features of the BP 3BU1-5**

The BP 3BU1-5 blood pressure monitor (with integrated time/date display) is a fully automatic, digital blood pressure measuring device for use on the wrist. This enables very fast and reliable measurement of the systolic and diastolic blood pressure as well as pulse by way of the oscillometric method of measuring. This device offers a very high and clinically tested measurement accuracy and is designed for maximum user-friendliness. Before using, please read through this instruction manual carefully and keep it in a safe place. For further questions on the subject of blood pressure and its measurement, please contact your doctor.

### **Attention!**

### **1.2. Important information about self-measurement**

Self-measurement means control, not diagnosis or treatment. Your values must always be discussed with your physician. You should never alter the dosages of any medication without direction from your physician. The pulse display is not suitable for checking the frequency of heart pacemakers!!

In cases of irregular heartbeat (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with your doctor.

Electromagnetic Interference: This device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens, etc.). These can lead to temporary impairment of the monitor's measuring accuracy.

### **2. Important information on the subject of blood pressure and its measurement**

### 2.1 How does high/low blood pressure arise?

Your level of blood pressure is determined in the circulatory center of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is effected by way of fine muscles in the blood vessel walls. Your level of arterial blood pressure changes periodically during heart activity. During the "blood ejection" (Systole), the value is highest (systolic blood pressure value). At the end of the heart's "rest period" (Diastole) pressure is lowest (diastolic blood pressure value). Your blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

### 2.2 Which values are normal?

Blood pressure is too high if your diastolic pressure is above 90 mmHg and/or your systolic blood pressure is over 160 mmHg, while at rest. In this case, please consult your physician immediately. Long-term values at this level endanger your health through continual damage to the blood vessels in you body. Should your systolic blood pressure values lie between 140 mmHg and 160mmHg and/or your diastolic blood pressure values lie between 90 mmHg and 95mmHg. Consult your physician. Regular self-checks will be necessary. If you have blood pressure values that are too low, (i.e., systolic values under 105mmHg and/or diastolic values under 60 mmHg), consult your physician. Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately. If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your physician. **Never use the results of your measurements to independently alter the drug doses prescribed by your physician.**  
**Which values are normal?**

Range	Systolic Blood Pressure	Diastolic Blood Pressure	Measures
Hypotension (low)	Lower than 100	Lower than 60	Consult your doctor
Normal Range	Between 100 and 140	Between 60 and 90	Self-check
Mild Hypertension (high)	Between 140 and 160	Between 90 and 100	Consult your doctor
Moderately Serious Hypertension (higher)	Between 160 and 180	Between 100 and 110	Consult your doctor
Serious Hypertension (very High)	Higher than 180	Higher than 110	Consult your doctor URGENTLY
Specific Systolic Hypertension (World Health Organization)	Higher than 140	Lower than 90	Consult your doctor

### Further information

- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called "labile hypertension." Consult your physician.
- Correctly measured diastolic blood pressure values above 120mmHg require immediate medical treatment.

### 2.3 What can be done if regular high/low values are obtained?

- Consult your physician.
- Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constrictions caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result. Furthermore, the heart will become structurally damaged with increased blood pressure values.
- There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. Secondary hypertension can be ascribed to specific organ malfunctions. Please consult your physician for information about the possible origins of your own increased blood pressure values.

d) There are measures you can take to reduce and even prevent high blood pressure. These measures must be permanent lifestyle changes.

1) Eating habits:

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

2) Previous illnesses

- Consistently follow all medical instructions for treating illness such as: diabetes (diabetes mellitus or sugar diabetes), fat metabolism disorder, and gout.

3) Habits


- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (coffee, tea and chocolate)

e) Your physical condition

- After a medical examination, and with your doctor's approval and direction; exercise.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your physician before beginning your exercise routine. You must receive advice regarding the type and extent of exercise that is appropriate for you.

### 3. The various components of your blood pressure monitor

The illustration shows your blood pressure monitor (model # BP 3BU1-5), consisting of:

Fits wrist size 13.5-19.5cm (5.3"-7.7")

### 4. Putting the blood pressure monitor into operation

#### 4.1. Inserting the batteries

After you have unpacked your device, insert the batteries. The battery compartment is located on the left lateral side of the device (see illustration).

- a) Remove cover as illustrated.
- b) Insert the batteries (2 x size AAA 1.5V), observing the indicated polarity.
- c) If the battery warning appears in the display, the batteries expired and must be replaced by new ones.

**Attention!**

- After the battery warning appears, the device is blocked until the batteries have been replaced.
- Please use «AAA» Long-Life or Alkaline 1.5V Batteries. The use of rechargeable batteries is not recommended.
- If the blood pressure monitor is left unused for long periods, please remove the batteries from the device.



**Functional check:**

Press and hold the O/I (ON/OFF) button to test all the display elements. When functioning correctly all segments appear.

**4.2. Sleep mode**

The device is designed to stay at the sleep mode as batteries are loaded. In this mode, the time and user icon set will be displayed. You can conduct the following operations during this mode.


- Reading the set date Press the TIME button. The date set will be displayed for 2 seconds.
- Recalling the measurements Press the MEMORY button. The measuring data stored can be displayed. For more information refer to the section on Memory storage and recall.
- Printing the measuring data Please press the PRINT button. The stored measurements can be printed. For details, please refer to the section on Printer Function.


**4.3. Setting the time and date**

This blood pressure monitor incorporates an integrated clock with date display. This has the advantage, that at each measurement procedure, not only the blood pressure values are stored, but also the exact moment of the measurement. After new batteries have been inserted, the clock begins to run from the following setting: 2001-01-01 00:00 O'clock. You must then re-enter the date and current time. For this, please proceed as follows (Example: Entering 2001-06-20 Time 09:30 o'clock):


- 1) Press and hold the Time button for three seconds. User 1 or 2 will begin flashing. Press the Time button again briefly. Now you can begin changing the time and date.


- 2) The correct year can be entered by pressing the MEMORY button. (Example: 1 x press)

- 3) Press the TIME button again. The display now switches to the current date, during which the first character (month) blinks.


- 4) The corresponding month can now be entered by pressing the MEMORY button. Example: 6 x presses)


5) Press the TIME button again. The last two characters (day) are now blinking


6) The corresponding day can now be entered by pressing the MEMORY button. (Example: 15 x presses)


7) Press the TIME button again. The display now switches to the current time, during which the first character (Hour) blinks


8) The corresponding hour can now be entered by pressing the MEMORY button. (Example: 9 x presses)


9) Press the TIME button again. The last two characters (Minutes) now blink

10) The exact time can now be entered by pressing the MEMORY button. (Example: 30 x presses)


11) Now after all settings have been made, press the TIME button once again. The date is briefly displayed and then the time. The input is now confirmed and the clock begins to run.


#### 4.4. User selection:

This advanced blood pressure monitor allows you to track blood pressure readings for 2 individuals independently.

a) Before measurement make sure you set the unit for the intended user. The unit can track results for 2 individuals. (User 1, User 2)


b) The unit is set to User 1. Push and hold the Time button for 4 seconds to initiate the user selection.


c) When the user icon is flashing, press the MEMORY button to toggle between users.

d) When the correct user is selected wait three seconds for the user icon to stop flashing. The monitor is now ready for use.

e) Repeat the procedure again to change users again.

#### Further Information

With each press of the button (TIME, MEMORY) one input is made (e.g. switching over from hours to minutes mode, or altering the value by +1). However, if you keep the respective button depressed, you can switch more quickly to the desired mode, or find the desired value respectively.

## 5. Carrying out a measurement

### 5.1. Before the measurement:

- Avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Try and find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before the measurement.
- Always measure on the same wrist (normally left).
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.

### 5.2. Common sources of error:


#### Note:

Comparable blood pressure measurements always require the same conditions! These are normally always quiet conditions.

- All efforts by the patient to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the wrist artery lies considerably lower (higher) than the heart, an erroneously higher (lower) blood pressure will be measured! (Each 15cm/6" difference in height results in a measurement error of 10mmHg!)
- A loose cuff causes false measurement values.
- With repeated measurements, blood accumulates in the respective arm, which can lead to false results. Correctly executed blood pressure measurements should therefore first be repeated after a 5 minute pause or after the arm has been held up in order to allow the accumulated blood to flow away (after at least 3 minutes).


### 5.3. Fitting the cuff

- a) Remove all objects and jewelry (e.g. wrist-watch) from the wrist in use. Put the cuff over the wrist with the monitor face on the inside of the wrist.

12

- b) The distance between the cuff and the hand should be approx. 10 mm (1/2").

- c) Secure the cuff with the Velcro fastener, so that it lies comfortably and not too tight. No space should remain between the cuff and the wrist.

- d) Lay your arm on a table, with the palm upwards.

Support the arm a little with a cushion, so that the cuff rests at about the same height as the heart. Remain still for 2 minutes, before beginning the measurement.

### 5.4. Measuring procedure

After the cuff has been appropriately positioned, and you've selected User 1 or User 2, the measurement can begin:

- a) If you want to take an Average Mode measurement position the slide switch to the 3. Remember, in this mode, the unit takes 3 measurements (one after another) and calculates the result. There is a 35 second rest time between measurements. The unit beeps five seconds before each upcoming measurement. (Data obtained in Average Mode are indicated by an A symbol.)

13




b) If you want to take a Single Mode measurement position the slide switch to the 1. In this mode there is only one measurement cycle.

c) Once you've selected the correct user, and measurement mode you can press the 0/I (ON/OFF) button to begin. After pressing the 0/I button the pump begins to inflate the cuff. In the display, the increasing cuff-pressure is continually displayed.

d) After reaching the inflation pressure, the pump stops and the pressure slowly falls. The cuff-pressure (large characters) is displayed during the measurement. When the device has detected the pulse, the heart symbol in the display begins to blink and a beep tone is audible for every pulse beat.


e) When the measurement has been concluded, a long beep tone sounds. The measured systolic and diastolic blood pressure values as well as the frequency are now displayed. (Pulse=number of heartbeats each minute) Example (Fig.): Systolic 118, Diastolic 73, Pulse 75. The measurement results are displayed, until you switch the device off. If no button is pressed for 5 minutes, the device switches automatically off, to save the batteries.

f) When the unit is set to Average Mode setting 3 separate measurements will take place in succession followed by a calculation of your detected blood pressure value. There will be 35 seconds of resting time in-between each measurement. A count down indicates the remaining time and a beep will sound 5 seconds before the 2nd and 3rd readings begin. If one of the measurements causes an error message, it will be repeated one more time. In case any additional error occurs, the measurement will be discontinued and an error code displayed.


#### 5.5. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the 0/I (ON/OFF) button can be pressed at any time. The device then immediately lowers the cuff-pressure automatically and enters into sleep mode.

### 5.6. Memory – storage and recall of the measurements

This blood pressure monitor automatically stores each of the last 30 measurement values per user. (Be sure to select the desired user (either user 1 or user 2) before accessing the stored measurements.) By pressing the MEMORY button, the last measurement (MR1) as well as the last 29 measurements of that user (MR2, MR3, ..., MR30) can be displayed one after the other. (MR1: Values of the last measurement.) (MR2-MR30: Values of the measurement before MR1.)


(MR1: Values of the last measurement.)

(MR2: Values of the measurement before MR1.)

#### Further Information

When the unit is in Average Mode, time between measurements is carefully controlled to ensure accuracy. If you take successive measurements when the unit is in the single measurement mode wait for several minutes in a relaxed position before you repeat a measurement. Make sure all measurements are taken in a seated position.

### 5.7. Memory-recall of 3 intermittent average results

The blood pressure monitor can display the 3 intermittent average results of the last stored memory. By pressing the MEMORY button at least 3 seconds, the LCD will show **A** icon first with one short beep.

Then release the MEMORY button, and the 3 intermittent average results will appear. NOTE: This only occurs if the last stored measurement is taken in Average Mode.

### 5.8. Memory– cancellation of all measurements


#### Attention!

To delete all stored readings, press the MEMORY button for at least 7 seconds. The display will show the symbol "CL" and 3 short beeps will sound to indicate deletion of stored readings.

#### Note:

If your last stored measurement was taken in the Average Mode, you will see **A** icon with a short beep during the deletion procedure. (The display will show these symbols during deletion of all stored measurements if the last measurement is taken from the Average Mode)

### 6. PRINTER FUNCTIONS

This unit can be used in connection with the Microlife diagnostic printer. The unit can print out measurements stored in the memory by pressing the PRINTER button. There are 2 kinds of formats available. Please make sure you set the unit to User 1 or User 2 based on the data you want to print.



**6.1. How to print the last stored measurement (Single data print)**

a) Select correct printer objective first before you press the PRINTER button.



b) Connect the printer with the monitor and switch the printer on (for details please refer to the manual of the printer).

c) Press the PRINTER button of the monitor once. The display will show "Pr1" and the printer will print the last stored measurement for User 1 or 2 based on your selection.


**6.2. How to print all measurements stored in memory and a data chart, (data print out with chart).**

a) Select correct printer objective first before you depress the PRINTER button.

b) Connect the printer with the monitor and switch the printer on (for details please refer to the manual of the printer).

c) Depress the PRINTER button for at least more than 3 seconds. The display will show "Pr1" followed by "Pr2"


d) Release the button, and the printer will print all stored measurements with a corresponding chart.

**6.3. Stop printing**

If you would like to stop printing, click the PRINTER button during the printing process.


**NOTE:**

Further details regarding the printer can be found in the printer instruction manual.

**7. Error messages/malfunctions**

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (Example: Error No. 2).

Error No.	Possible cause(s)
ERR 1	No pulse has been detected.
ERR 2	Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the Measurement (Artefact).
ERR 3	The inflation of the cuff takes too long. The cuff is not correctly seated.
ERR 5	The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take another reading following directions carefully. Contact your doctor if you continue to get unusual readings.

### Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains blank when the instrument is switched on although the necessary insert correctly.	1. Check batteries for correct polarity and if batteries are in place. 2. If the display is unusual, re-insert batteries or exchange them.
The device frequently fails to measure the blood pressure values, or the values measured are too low (too high).	1. Check the positioning of the cuff. 2. Measure the blood pressure again in peace and quiet under observance of the details made in Section 5.
Every measurement produces a different value	1. Please read the following information and listed in Section 5.2 " <b>Common sources of error.</b> " Repeat the measurement. <b>Please note: Blood pressure fluctuates continually so successive measurements will show some variability.</b>
Blood pressure values measured differ from those values	1. Record the daily development of the measured by the doctor and consult your doctor. <b>Please note: Individuals visiting their doctor frequently experience anxiety which can result in a higher reading at the doctor than obtained at home under resting conditions.</b>


### Further Information

Your level of blood pressure is subject to fluctuations even with healthy people. Comparable measurements always require the same conditions (Quiet conditions)! If, in spite of observing all these factors, the fluctuations are larger than 15mmHg, and/or you hear irregular pulse tones, please consult your doctor. For licensing, this device has been subjected to strict clinical tests. The computer program used to measure the blood pressure values was tested by experienced specialist physicians in Germany. The same computer program is used in every individual device, and has also been clinically tested. The manu-

facture of the devices takes place according to the terms of the European standard for blood pressure measuring devices (see technical data) under the supervision of the Technical Monitoring Association Essen (RWTÜV-Essen). You must consult Microlife if there are technical problems with the blood pressure instrument. Never attempt to repair the instrument yourself! Any unauthorised opening of the instrument invalidates all warranty claims!

### 8. Care and maintenance, recalibration

- Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.
- Clean the device with a soft, dry cloth. Do not use gasoline, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. The cuff must not be washed!
- Do not drop the instrument or treat it roughly in anyway. Avoid strong vibrations.
- Never open the device! Otherwise the manufacturer calibration becomes invalid!

### Periodical recalibration

Sensitive measuring devices must from time to time be checked for accuracy. We therefore recommend a periodical inspection of the static pressure display every 2 years.

Microlife is pleased to provide more extensive information on calibration.

### 9. Warranty

The blood pressure monitor BP 3BU1-5 is guaranteed for 5 years from date of purchase. This warranty includes the instrument and the cuff. The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or alterations made to the instrument by third parties. The warranty is only valid upon presentation of the warranty card.

### 10. Reference to standards

#### Device standard:

Device corresponds to the requirements of the European standard for non-invasive blood- pressure monitor  
 EN1060-1 / 12:95  
 EN1060-3 / 09:97  
 DIN 58130, NIBP – clinical investigation  
 ANSI / AAMI SP10, NIBP – requirements  
 (US Standard)

#### Electromagnetic compatibility:

Device fulfills the stipulations of the European standard EN 60601-1-2

#### Clinical testing:

The clinical performance test was carried out in Germany according to the DIN 58130 / 1997 procedure N6 (sequential). The stipulations of the EU-Guidelines 93/42/EWG for Medical Products Class IIa have been fulfilled.

### 11. Technical specifications:

Weight:	148g (with batteries + cuff)
Size:	85 x 77 x 75mm (including cuff)
Storage temperature:	-5 to +50°C (23°F - 122°F)
Humidity:	15 to 85% relative humidity maximum
Operation temperature:	10 to 40°C (50°F - 104°F)
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	oscillometric
Pressure sensor:	capacitive
Measuring range:	SYS/DIA: 30 to 280 mmHg
Pulse:	40 to 200 per minute
Cuff pressure display range:	0-299 mmHg
Memory:	Stores the last 30 measurements automatically for 2 Users
Measuring resolution:	1 mmHg
Accuracy:	Pressure within ± 3 mmHg Pulse ± 5 % of the reading
Power source:	2 dry cells (Batteries) UM-4, size AAA 1.5V
Accessories:	storage case
Technical alterations reserved!	

### 12. [www.microlifeusa.com](http://www.microlifeusa.com)

Detailed user information about our products as well as services can be found at [www.microlifeusa.com](http://www.microlifeusa.com) or call us toll free at 800-568-4147. Microlife USA, Inc.  
 26Morin St  
 PO Box 530  
 Biddeford, ME 04005  
 Telephone: 207-284-1012  
 Fax: 207-283-1122  
 Email: [custserv@microlifeusa.com](mailto:custserv@microlifeusa.com)

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