microlife Thermal Printer

Model# PR 1KA1 For Premier Automatic Blood Pressure Monitors
(#3AA1-2 and #3BU1-5) and Peak Flow Meter

Instruction Manual

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microlife[®]

PLEASE NOTE:
THIS MEDICAL INSTRUMENT MUST BE USED ACCORDING TO
INSTRUCTIONS TO ENSURE ACCURATE READINGS.

Questions: Call microlife toll-free at 1-800-568-4147

Diagnostic Printer

Instruction Manual

Table of contents

1. Introduction

2. Printer parts

3. Putting the printer into operation

- 3.1. Inserting the batteries
- 3.2. Using an AC adapter
- 3.3. Cable connection
- 3.4. Inserting the paper roll

4. Carrying out a printing operation

- 4.1. Power on
- 4.2. Paper feeding
- 4.3. Printing
- 4.4. Low voltage indication
- 4.5. Power off
- 4.6. Cutting paper

5. Malfunction and Remedy

6. Care and maintenance

7. Warranty

8. International quality standards

9. Technical specifications

10. www.microlife.com


1. INTRODUCTION

The diagnostic printer PR 1KA1 is a compact and small unit for Microlife Premier Blood Pressure Monitors and Asthma Monitors. The enhanced printer covers both character and graphic print. High-reliability and high-quality print is assured due to the ultra simple design and practicability. Key features are:


- Compact and lightweight
- High-quality
- Battery-powered
- Designed for versatile applications
- Simple mechanism ensures high-reliability

2. PRINTER PARTS

a) Front View

b) Back View


3. PUTTING THE PRINTER INTO OPERATION

3.1. Inserting the batteries


The battery compartment is located on the underside of the device (Please see illustration).

- Remove cover as illustrated.
- Insert the batteries (4 x size AA), thereby observing the indicated polarity.
- If the LED starts flashing, the batteries are empty and must be replaced.

Attention!

- If the green LED starts flashing, the device is blocked until the batteries have been replaced.
- Please use 'AA' Long-Life or Alkaline 1.5V Batteries. Do not use rechargeable batteries.
- If the printer is left unused for long periods, please remove the batteries from the device.


3.2. Using an AC adapter

It is possible to operate this printer through the Microlife AC adapter (output 6V DC / 600mA with DIN plug). Make certain that you use an AC adapter which fulfills the legal requirements.

(CE sign on the type plate).

- Plug the DIN plug into the socket at the left back of the instrument.
- Plug the AC adapter into a 110V power socket. Test that power is available by sliding the O/I (on/off) switch.

Note:


- No power is taken from the batteries while the AC adapter is connected to the instrument.
- If the power is interrupted during the measurement (e.g. by removal of the AC adapter from the socket), the instrument must be reset by removing the DIN plug from its socket.
- Please consult Microlife if you have questions relating to the AC adapter.

3.3. Cable connection

Insert one end of the cable into the printer port socket (right backside of the instrument) as shown in the diagram.




Note:

- Connect the other end with the monitor or asthma monitor (both connectors are identical).

3.4. Inserting the paper roll


- Switch the printer on.
- Insert the end of paper roll into the slot of the internal compartment.
- Move the paper as deep as possible into the slot.
- Press the paper feed button until paper appears at the front.
- Put paper roll into the compartment and close the lid (small paper roll) or put the paper roll into the lid (big paper roll).

There are 2 kinds of paper rolls for use. The small size of paper roll is placed into internal compartment. The big size of paper roll has to be placed into the outer compartment which is the inside of the printer cover. Microlife sells extra paper rolls.

Note:

- The small paper roll is 1" thick (25mm-standard package) and 38mm wide.
- The big paper roll is 2" thick (50mm) and 38mm wide.
- If there are problems when inserting the paper, cut the edge of the paper so that the corners of the paper are square (90°, see below).

4. CARRYING OUT A PRINTING OPERATION

4.1. Power on

Before powering the printer, check if batteries are inserted properly.

- Pushing the slide switch to "left" side in order to power on the printer.
- Simultaneously, the green LED will become bright.


4.2. Paper feeding

Before you start the printing operation, check if the paper is installed properly in the right position. If it is not, click the paper feed button to get access to proper location.

Note:

- Click once to get 1 line space.
- Continual depressing to get multi lines space.

4.3. Printing

- Connect the printer with the monitor as explained at 3.3.
- Push the print button on the monitor (for details please refer to the manual of related monitor).

4.4. Low voltage indication

Once the LED is flashing under power-on operation, it indicates that new batteries are needed.

- Once the batteries are replaced the printer will work.

Note:

- If the printer is connected with an adapter, this will not happen.

4.5. Power off

After the printing task, just push the slide switch to the "right" side in order to power off the printer.

- Simultaneously, the LED will become dark.


4.6. Cutting paper

Once you finish the printing job, check the paper to see if it is easy to tear off.

- If it is not long enough to tear off, please depress the "Paper Feeding" button to lengthen the paper.

Note:

- Don't pull out the paper if you'd like to cut it.
- Try to lower down the paper, pulling it forward over the paper cutter.
- The paper has to be vertical (90 degree) with the paper cutter; otherwise, the paper will be difficult to tear off.
- The print on the paper does not last forever. It fades over time. This is normal for thermal printers such as the Microlife diagnostic printer. The readability lasts longer the less the print is in contact with light. If you file the printed paper in a folder, the print lasts very long. In case you would like to store the values forever, we recommend to photocopy them or transfer the data to another media.

5. MALFUNCTION AND REMEDY

If problems occur when using the printer, the following points should be checked and if necessary, the corresponding remedies taken:

Malfunction	Remedy
The LED remains flashing while the printer is switched on although the batteries are good.	<ol style="list-style-type: none">1. Check whether the batteries are installed with the correct polarity.2. If the LED is unusually flashing, remove the batteries and re-install them or exchange them for new.
Clicking the printing button on the monitor, but nothing prints	<ol style="list-style-type: none">1. Check whether the cable is connected with the monitor and printer.2. Examine whether the paper is positioned properly ready for printing.
There's no response for paper feeding	<ol style="list-style-type: none">1. Examine whether the paper is positioned properly.2. Check if the paper roll has been exhausted, and needs to be replaced by a new one.

6. CARE AND MAINTENANCE

- a) Do not expose the device to extreme temperatures, humidity, dust or direct sunlight.

- b) Clean the device with a soft, dry cloth. Do not use gasoline, thinners or similar solvent.

- c) Do not drop the device or treat it roughly in any way. Avoid strong vibrations.

- d) **Never open the device!**
Otherwise the manufacturer warranty is invalid.

7. WARRANTY

The diagnostic printer PR 1KA1 is guaranteed for **2 years** from date of purchase. The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or alterations made to the device by third parties. The warranty is only valid upon presentation of the completed warranty card. Batteries or damage from batteries are not covered under warranty.

Name and company address of the responsible dealer:

8. INTERNATIONAL QUALITY STANDARDS

Meets the requirements of MDD93/92/EEC (Europe) and U.L. (USA).

9. TECHNICAL SPECIFICATIONS:

Weight:	255g (with batteries)
Size:	143 (L) x 81 (W) x 51 (H) mm
Storage temperature / Humidity:	-20 to +50°C / 15 to 90% relative humidity maximum
Operation temperature / Humidity:	10 to 40°C / 15 to 90% relative humidity maximum
Print system:	Thermal serial dot system
Print direction:	Left to right
Printer head model no.:	SII MTP102-16B
Printer head dimension:	48(L) x 31(W x 13.8(L)
Printing timing:	Synchronized with the tachogenerator
Home position detection:	Mechanical switching system
Printing speed:	Approx. 1.5 (lines/sec.)
Number of columns:	26 (columns)
Character matrix (HxW):	7x5 (dot)
Printing width / Paper width:	29.0 / 38.0 (mm)
Reliability:	500,000 (lines)
Paper roll:	1" 25mm (standard size); 2" 50mm (large size)
Power source:	1.5V (AA) alkaline batteries x 4
Low battery:	As battery is < 5.0V, LED will flash continuously with 3 Hz
Accessories:	a) 100cm cable b) Instruction manual
Optional:	AC/DC adapter 6V DC, 600mA, 50-60 Hz Batteries sold separately

* Technical alterations reserved!

10. www.microlife.com

Detailed user information about our products as well as services can be found at www.microlife.com

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Microlife USA, Inc.