60-Second Basal Thermometer

INDICATIONS FOR USE:

This clinical digital thermometer is a device intended for measuring, and aiding in the monitoring and tracking of basal body temperature (BBT) as an aid in ovulation prediction for conception. Tracking your BBT helps predict ovulation or fertility timing by predicting maximum and minimum fertile days to help increase the chances of conception or pregnancy.

This thermometer IS intended to be used:

- For measuring, and aiding in the monitoring and tracking of basal body temperature (BBT) as an aid in ovulation prediction for conception.
- Within a home healthcare environment
- Oral & Vaginal
- With an optional Probe Cover (probe cover not mandatory)

This thermometer IS NOT intended for use with:

- Contraception

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE USE

- Do not use this device if you think it is damaged or anything appears unusual
- To avoid possible injury, do not run or walk while using this device.
- Ensure that children are supervised around this medical device. Keep batteries away from children; harmful if swallowed.
- Do not dispose of the battery or device in fire; the battery may explode. Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations as domestic waste.
- Do not use a metal tool to replace the battery; doing so may short circuit the terminals.
- Do not attempt to recharge the battery as they are not rechargeable.
- Only open the thermometer battery cover to replace the battery (when applicable); do not open/modify this device; inaccuracy and/or harmful injury may result.
- Performance of this thermometer may be degraded if used outside of its specified operation temperature and humidity rating, or if stored outside of its specified storage temperature and humidity rating (avoid storage in direct sunlight), the device is exposed to mechanical shock, the patient's temperature is lower than that of the ambient, the patient's temperature is outside the specified patient range, the tip is dirty or damaged (ensure the tip is clean before using), the optional probe cover is dirty or a probe cover not intended for use with this device is used.
- Do not expose the thermometer to extreme temperature, humidity, shock, and/or dust.
- Any probe covers that may be included in this device package are for single use only.
- This device contains sensitive electronics components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g., mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring accuracy. Move the device to another location if interference is determined and maintain a minimum distance of 11 ft from such devices when using this unit.
- It is recommended that this digital thermometer be checked for calibration every 2 years

WHAT IS A BASAL TEMPERATURE?

A basal temperature is the normal body temperature of a healthy person immediately upon awakening after a restful night's sleep. More precisely, it is the body temperature measured under so-called basal conditions (12 hours after eating, after a restful sleep, no exercise, no emotional excitement, normal room temperature).


For optimum accuracy, your basal temperature should be taken at the same time every morning upon waking. Changes in basal temperature are related to changes in the female reproductive cycle. By taking your basal temperature daily and charting it, you can notice slight changes and patterns. These charts of your basal temperature, in addition to other information, can help you understand your menstrual cycle and ovulation. Recording and interpreting basal temperature changes is called the Basal Temperature Method. This method can be used to help in family planning. Charting of temperature and prediction of ovulation can be used to time sexual intercourse during fertile days to aid in increasing the likelihood of pregnancy.

HOW TO TAKE A TEMPERATURE:

This thermometer can be used with or without a probe cover. To use a cover, slide onto sensor tip until fully engaged. If thermometer has been in a cold environment, allow it to warm up to room temperature for about 15 minutes before use.


Avoid storing at temperatures above 140°F (60°C).

To turn the unit on, press the power button. Proper functioning of the thermometer is tested automatically each time it is turned on. The display should read "188.88 °E" to confirm normal function. An "Err" message indicates unit is not working properly. The display will then show "Lo" and a flashing "F" or "C" indicating that thermometer is ready for use.


ORAL USE:

- Do not drink hot or cold liquids or exercise before use.
- Place sensor tip under tongue as shown in Figure A.
- Accurate temperature reading will take approximately 60 seconds. An open mouth or improper placement could make reading take longer.

VAGINAL USE:

Please seek guidance from your health care professional. Taking vaginal temperatures is beyond the scope of these instructions.

READING THE TEMPERATURE:

Hormonal changes associated with ovulation produce a slight rise in a woman's basal body temperature. A change of approximately 0.5°F (0.3°C) occurs in approximately the last two weeks of the cycle around the time of ovulation. This temperature shift, as well as other recognizable signs such as changes in cervical mucus, can assist in determining the fertile days prior to ovulation.

CHANGING FROM °F TO °C:

With the thermometer turned off, press and hold the power button approximately 2 seconds until °F or °C changes on the display. The thermometer will briefly display a self test temperature. It will then be set to your selection, ready for use.

HOW TO CLEAN:

Clean thermometer after every use:

1. Clean sensor tip with 70% isopropyl alcohol. Wipe dry with sterile cotton ball.

HOW TO CHANGE BATTERY:

Replace battery when "▽" appears on display.

1. Remove the battery cover by pulling away from the base.
2. Using a toothpick (or other non-metal object) remove the old battery.
3. Insert new 1.5 V button cell, LR41 or equivalent, battery with the contact down (+ symbol facing toward you).
4. Replace battery cover.

LIMITED WARRANTY:

This thermometer is warrantied for one year by Microlife USA Inc, against manufacturer defects for the original purchaser only, from date of purchase. The one year warranty applies to the thermometer. Batteries and case are not covered by this warranty. There are no user serviceable parts inside.

The warranty does not apply to consequential and incidental damages, or damage caused by batteries, improper handling, and accidents. Professional use, not following the operating instructions, and alterations made to the device or accessory by third parties, are also not included in this warranty. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

Microlife USA Inc will investigate your concern. A product as defined by this warranty, determined to be out of specification, will be replaced and shipped to you at no cost. A product or accessory as defined by this warranty, determined to be within specification, will be returned to you with a report of findings, at no cost. Please call customer service: 1-800-568-4147 regarding any warranty concerns. We ask that you please contact us before sending any product back in order to better identify, and more quickly process, your concern.

TECHNICAL SPECIFICATIONS:

Measurement range:	90.00°F (32.00°C) to 109.99°F (42.99°C)
Measurement accuracy:	±0.2°F (±0.1°C) between 90.00°F (32.00°C) to 109.99°F (42.99°C)
Self test:	Automatic internal check of sensor circuit. If error, "Err" is displayed.
Display:	Liquid crystal display (LCD) with five digits
Auto Off:	After 10 minutes of inactivity
Smallest unit of display:	0.01°F (0.01°C)
Memory:	Stores the last measured value
Operating conditions:	50°F (10°C) to 104°F (40°C), humidity of 15% to 95% (non-condensing)
Air pressure:	700 to 1040 hPa
Storage conditions:	-13°F (-25°C) to 140°F (60°C), humidity of 15% to 95% (non-condensing)
Battery:	1.5V, LR41
Product life:	5 years
Reference to standards:	Device corresponds to the requirements of the standard for electronic thermometers ASTM E-1112, EN 12470-3, IEC 60601-1, 60601-1-2, 60601-1-11. Complies with ISO 10993-10:2010, ISO 10993-5: 2009, ISO 10993-12:2012
Device standard:	
IP22:	Protected against solid foreign objects of 12.5 mm diameter and greater. Protected against vertically falling water drops when the device is tilted up to 15°.
Electromagnetic Compatibility:	Device fulfills the stipulations of the standard IEC 60601-1-2

Test Protocol (for laboratory use only):

1. Place thermometer into waterbath. The tip of the thermometer should be approximately 2cm below the surface.
2. Wait for 10 seconds, then turn the thermometer on by pressing the "On/Off" button. The thermometer is now in standard mode.
3. When the measurement is completed, record the temperature reading.

Distributed by:
Microlife USA, Inc.
1617 Gulf to Bay Blvd., 2nd Floor
Clearwater, FL 33755
©2020 Microlife Corporation
www.microlifeusa.com

Questions?
We're here to help. Call toll-free
1-800-568-4147

Follow instructions of use to avoid hazards.

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