# microlife®

## INSTRUCTION MANUAL

Digital  
Peak Flow &  
FEV1 Meter

Model: PF100

**QUESTIONS? 1-866-934-1839**

# Digital Peak Flow Meter – Instruction Manual

## Table of Contents

---

### **1. Introduction**

---

### **2. Important Information Concerning Intended Use**

- 2.1 Information about the Measurement and Control of Values
- 2.2 Which values are Normal?

---

### **3. Components of your Peak Flow Monitor**

---

### **4. Using your Monitor for the First Time**

- 4.1 Inserting the Batteries
- 4.2 Setting Date and Time

---

### **5. Measurement procedure**

- 5.1 Measuring Procedure
- 5.2 Information about the Measurement and Control of Peak Flow Values
- 5.3 Self-Assessment with the Traffic Light Indicator
- 5.4 Stored Readings/Memory
- 5.5 Battery Replacement

---

### **6. Software**

---

### **7. Error Messages/Troubleshooting**

---

### **8. Care and Maintenance**

- 8.1 Cleaning / Disinfecting the Mouthpiece & Measuring Tube
- 8.2 Cleaning the Main Unit
- 8.3 Maintenance and Safety

---

### **9. Limited Warranty**

---

### **10. Certifications**

---

### **11. Technical Specifications**

---

### **12. How to Contact Us**


Dear Customer:

Please use this Traffic Light card following your physician's (or licensed health care professional's) advice. They will help determine good flow rates and recommend actions when decreased flow rates are measured.

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## **1. Introduction**

---

Your new Microlife Peak Flow Meter is a high quality medical device that measures your maximum possible exhalation which is called peak flow (referred to in medical terms as PEF – Peak Expiratory Flow) and forced expiratory volume (commonly known as FEV1) in 1-second. Regular monitoring of your peak flow is extremely useful for controlling diseases of the airways such as asthma, C.O.P.D., or chronic bronchitis.\*

The Microlife monitor is very simple to use, so measurements can easily be performed at home. The Peak Flow Meter is designed for people of all ages and is suitable for all ages. Additionally, it automatically stores 240 readings so you can take your daily results with you when you visit your doctor for a checkup.

Read through this instruction manual carefully before using the device and then keep it in a safe place.

*\*A higher value usually means air is moving easily through your lungs. When asthma attacks occur air cannot move easily through your lungs and lower values will be noted.*

## **2. Important Information Concerning Intended Use**

---

### **2.1 Information about the Measurement and Control of Peak Flow Values**

This device is intended for monitoring PEF (Peak Expiratory Flow Rate) and FEV1 (Forced Expiratory Volume in one second) for patient home use. The device is designed for pediatric to adult patients, and the simple interface makes it easy to use.

When the device is used to monitor lung conditions such as asthma and chronic obstructive pulmonary disease (COPD), the user should be under the care of a licensed health care

professional. A licensed health care professional's advice is required to understand the meaning and importance of the measurements reported by the device and how to decide on an appropriate treatment plan that defines when to measure.

A peak flow meter is used to measure a person's "peak expiratory flow," which is the fastest speed a person can blow air out of the lungs after taking in as big a breath as possible. "Peak expiratory flow" is a simple measure of airflow that can tell you how well you are breathing. It tells you how well air is moving through the airways in your lungs.

Forced Expiratory Volume (FEV1) is a measure of the volume of air expelled in 1-second.

If you have a breathing condition such as asthma or C.O.P.D., your physician (or other licensed health care professional) may recommend that you use a peak flow meter to watch for changes in your airflow.

Your physician (or licensed health care professional) will give you a treatment plan that will tell you what actions to take when you have a change in airflow. In addition, you should record your peak flow measures (see sections 5.4, Stored Readings/Memory) as recommended by your physician (or other licensed health care professional). Reviewing peak flow measures can help you and your physician (or licensed health care professional) check closely on your asthma or C.O.P.D. to provide the best treatment for you.

The treatment plan given to you by your physician or other licensed health care professional will tell you what action to

take when there are changes in your peak flow number.

No matter what your peak flow measures are, if you have signs and symptoms such as chest tightness, shortness of breath, coughing or wheezing you should follow your licensed health care professional's advice for contacting him or her.

Instructions for using peak flow meter must be followed carefully in order to get a correct measure of airflow.

If you are unable to obtain a reading, and your monitor is working correctly, you should contact your physician immediately.

## 2.2. Which PEF Values are Normal?

**Note:** Please contact your doctor for FEV1 standards.

Normal PEF Values for Males\*

Age (Years)	Height				
	60"	65"	70"	75"	80"
20	554	602	649	693	740
25	543	590	636	679	725
30	532	577	622	664	710
35	521	565	609	651	695
40	509	552	596	636	680
45	498	540	583	622	665
50	486	527	569	607	649
55	475	515	556	593	634
60	463	502	542	578	618
65	452	490	529	564	603
70	440	477	515	550	587

## Normal PEF Values for Females\*

Age (Years)	Height				
	55"	60"	65"	70"	75"
20	390	423	460	496	529
25	385	418	454	490	523
30	380	413	448	483	516
35	375	408	442	476	509
40	370	402	436	470	502
45	365	397	430	464	495
50	360	391	424	457	488
55	355	386	418	451	482
60	350	380	412	445	475
65	345	375	406	439	468
70	340	369	400	432	461

## Normal PEF Values for Children and Adolescents\*


Height (inches)	Males & Females
43"	147
44"	160
45"	173
46"	187
47"	200
48"	214
49"	227
50"	240
51"	254
52"	267
53"	280
54"	293
55"	307
56"	320
57"	334
58"	347
59"	360
60"	373
61"	387
62"	400
63"	413
64"	427
65"	440
66"	454

\*Leiner GC. et al: Expiratory peak flow rate. AM Rev Respir Dis 88:644, 1963

1 Polgar G. Promadhat V: Pulmonary Function Testing in Children: Techniques and Standards. Philadelphia, W.B. Saunders Company, 1971

### **3. Your Microlife Peak Flow Meter**

---

## **4. Preparation for your first Measurement**

---

Before your first use we recommend that you disinfect the mouthpiece and flow tube as explained in this manual. (Section 8)

### **4.1. Inserting the Batteries**


1. Please open the Battery Compartment from underneath and insert the two batteries (1.5 V, size AAA).
2. Follow the polarity as indicated by the symbols in the compartment.

## 4.2. Setting Date and Time

After new batteries have been inserted, the default time/date display will blink. You must then re-enter the current date and time. To do this, proceed as listed below. Once set, the monitor automatically records the time and date of each measurement:

1. Lift and open the data compartment from the side of the device as indicated by the data symbol.
2. If the year is not blinking in the display, press the small TIME button next to the clock symbol with a paper clip and the year will start blinking; release the switch.
3. Using the two arrow keys on the front of the monitor you can decrease (left arrow button) or increase (right arrow button) the number. Release the arrow button when the correct number is reached.
4. To confirm the year and then set the month, press the TIME button.
5. You can now set the month using the arrow buttons. (Example: pressing 2 times the right arrow button advances to 06 for June). To confirm and then set the day, press the TIME button.

6. Please follow the previous sequence to set the day, hour and minutes.
7. Once you have set the last minute and pressed the TIME button, the date and time are set and the time is displayed.
8. Close the data compartment.
9. If you want to change the date and time, press the TIME button and repeat steps 3-8.

**Please note:**

- If you hold down the arrow key for more than 2 seconds it speeds up.
- Date /time can also be set easily from computer when you run the Microlife Analyzer Software.

## **5. Measurement procedure**

---

**CAUTION: ASK YOUR PHYSICIAN (OR OTHER LICENSED HEALTH CARE PROFESSIONAL) TO WATCH YOU USE THIS PEAK FLOW METER BEFORE YOU RELY ON ANY MEASUREMENTS.**

**READ ALL INFORMATION IN THIS INSTRUCTION BOOKLET BEFORE USING.**

**CONTACT MICROLIFE AT 866-934-1839 FROM 8 AM – 5 PM (EASTERN TIME) MONDAY – FRIDAY WITH QUESTIONS.**


**IF YOU HAVE SIGNS AND SYMPTOMS SUCH AS CHEST TIGHTNESS, SHORTNESS OF BREATH, COUGHING OR WHEEZING FOLLOW YOUR PHYSICIAN'S (OR OTHER LICENSED HEALTH CARE PROFESSIONAL'S) ADVICE FOR CONTACTING THEM NO MATTER WHAT READINGS YOU OBTAIN.**

**IF YOU DO NOT GET A READING, AND YOUR DEVICE IS SET CORRECTLY (DATE AND TIME ARE SET, TUBE CORRECTLY IN PLACE) CONTACT YOUR PHYSICIAN IMMEDIATELY.**

**THIS MONITOR IS DESIGNED FOR A SINGLE USER.**

## 5.1. Measuring Procedure

1. Press the ON/OFF button to turn the monitor on. The last result will appear on the display (0 if there is no data). The device will sound two short beeps and arrows will blink to indicate it is ready for measurement.
2. You can perform the measurement while standing or sitting upright. For better comparison of your data, you should always perform the measurement in a similar position.
3. Hold the monitor with both hands on the two textured compartments.
4. Inhale completely and hold your breath for a moment.
5. Cover the mouthpiece tightly with your lips.
6. Blow into the measuring tube as hard and as fast as you can.
7. A long beep confirms that the result is recorded. PEF is displayed for about 3 seconds followed by FEV1. Then, the monitor will sound two short beeps indicating that it is ready for a new measurement.

8. It is recommended to perform three or more measurements sequentially.
9. After each measurement, the actual reading is shown followed by the highest reading of your current measurement session.
10. Press the ON/OFF button again to switch the monitor off. Before the monitor switches off, the highest value of the current measurement session and its related memory position MR XX is displayed.
11. Clean the tube and mouthpiece after use (see section 8).

**Please note:**

- **The monitor only saves the highest reading from a measurement session** (--measurements taken in one on/off cycle).
- The date and time must be set for the unit to take a measurement.
- The device stores up to 240 measured values with date and time. **When the memory is full, the oldest values are automatically erased. Make sure you visit your doctor in time to analyze your data.**
- Please confirm that date and time settings are always actual.
- The performance of the monitor can be affected by extremes of temperature and humidity. See technical specifications section for details.

## **5.2. Information about the Measurement and Control of Peak Flow Values**

- It is necessary to keep regular control of your peak flow values. Doctors generally recommend taking one measurement in the morning and one in the evening, every day, normally before inhaling your prescribed drugs. Carry out additional measurements whenever you feel unwell or perceive shortness of breath.
- Your doctor is interested in the highest peak flow value that you can achieve during a measurement session. Please repeat the measurements until you feel you have got the best possible result for the moment. Take care to have sufficient relaxation between measurements. In case of worsening results at each successive measurement session talk to your doctor. It may be a sign of unstable asthma.
- **Please note that asthma attacks are indicated in advance by low peak flow values before you feel it! Please contact your doctor in cases such as these. Also contact your doctor if you have signs and symptoms such as chest tightness, shortness of breath coughing or wheezing.**
- **Self-measurement is a means for control, and is not intended for diagnosis or treatment. In any event, please be sure to discuss your measured values with your doctor. Your doctor will also explain which values are normal for you.**
- You should never alter the dosages of any medication without talking to your doctor

### **5.3. Self-Assessment with the Traffic Light Indicator**

---

The colored bars on the left hand edge of the display are a quick visual reference to show the range of your measured values. As your peak flow raises or lowers, the arrow on the display will show the measurement as either within the normal (green), borderline (yellow) or danger (red) range.

As you use the monitor, it will automatically record your best value and use it to determine your traffic light zones.

#### **Determining the Area Limits**

The traffic light indicator can be used in one of two ways:

- Default traffic light zones
- Customized traffic light zones.


#### **Default Traffic Light Zones:**

1. The green zone ranges between your personal best measurement (MAX) and 80% of that reading.  
To determine the yellow/green border zone, multiply the MAX value by 0.8 [Example:  $500 \text{ l/min} \times 0.8 = 400 \text{ l/min}$ ].
2. The yellow zone ranges between your personal best measurement (MAX) and 60% of that reading.  
To determine the yellow/red border zone, multiply the MAX value by 0.6. [Example:  $500 \text{ l/min} \times 0.6 = 300 \text{ l/min}$ ].
3. The red zone ranges below 50% of your best reading.

#### **Customized Traffic Light Zones:**

The Traffic Light Indicator zones can be customized with a computer and the Asthma software. This Windows compatible software is free for download from our website (See section 6).

**CAUTION:** Setting area limits, creating a treatment plan and modifying your traffic light zone settings should only be done under direction of a physician (or other licensed health care practitioner). If you would like to use this feature, determine your customized zones with your doctor and record them on the traffic light card (enclosed).

**Note:** Once customized zones are programmed, default traffic light measurements will no longer be taken.

### Traffic Light Indicator Colors

The traffic light arrow will raise or lower according to the zones set for the monitor (default or customized).

#### Green Zone – OK

- If your measurement is in the green zone, your condition appears to be under control. Continue your treatment plan.

#### Yellow Zone – Caution


- If your measurement is in the yellow zone, measure more frequently and follow your treatment plan.

#### Red Zone – Danger

- If your measurement is in the red zone, your condition is serious. Act as discussed with your physician or seek emergency medical treatment.

## **5.4. Stored Readings/Memory**

Multiple readings can be taken during a session, and the monitor will automatically store the highest reading from the session (1 session=1 on-off cycle). This unit stores 240 readings.

### **Viewing the stored values**

1. Switch the monitor on and press the memory scroll buttons.
2. By pressing the < button, you can see the latest value stored in the memory with date and time. Continue to depress the button to move quickly through the memory.
3. The > button works in the opposite direction.

### **Memory Capacity Low**

When the memory reaches 230 or more readings, MR XX will blink after the monitor is switched on to indicate that the remaining memory capacity is low.

### **Memory Full**

When the memory has stored 240 results, the monitor will beep when the unit is switched on and MR 240 will blink. From this point, a new measured value will be stored by overwriting the oldest value.

## **Clear all Values**

To clear all data press the < and the > buttons simultaneously for 5 seconds, and release the buttons. CLR will appear and blink in the display for 3 seconds. To erase all data from memory, press the ON/OFF button while CLR is blinking. Otherwise the monitor moves back to normal operation and the data remains.

## **5.5. Battery Replacement**

---

When the battery low symbol appears in the display, the device is blocked until the batteries have been replaced

Use two 1.5 V Long-Life AAA batteries. If the monitor is left unused for longer periods, remove the batteries.

### **Notes:**

- During a battery exchange the data stored in the memory of the monitor is saved and will NOT be lost.
- The use of rechargeable batteries is NOT recommended.
- **After a battery change the time/date need to be reentered so that data can be stored.**

## **6. Software**

---

This meter can be used in connection with your personal computer (PC) by connecting the monitor via the included USB cable.

The free downloadable Microlife Asthma Software will enable you to transfer your monitor's readings to your PC, and to customize your traffic light zones (if desired). Once your readings are transferred, you will be able to print your results, view graphs and reports of your values and share them with your doctor.

**System Requirements:** Before installing your Asthma Monitoring Software, be sure that your computer meets or exceeds the system requirements:

<b>Hardware</b>	<b>Minimum</b>	<b>Recommended</b>
Microprocessor (CPU speed)	550 MHz	1 GHz or higher
Memory devices (RAM)	256 MB	512 MB
Hard disk	500 MB	800 MB
Communications port	USB 1.0	USB 1.0 and above
Display	256 colors	65536 colors
Resolution	800 x 600	1024 x 760
Operating system: Microsoft Windows XP, Vista, 7 & 8 Not compatible with macintosh computers.		

## **6.1. Installation and data transmission**

- a. Connect the monitor via USB cable with the PC. Please pay close attention when connecting the USB cable. Do not force the cable into the monitor or your computer.
- b) Windows OS – Install and launch the software. Please call customer service at 1-866-934-1839 for assistance.  
Microlife's Asthma Monitoring Software is available for download free of charge from our website:

**<https://www.microlife.com/consumer-products/respiratory-care/asthma-monitor/pf-100/>**

- c) Turn the Peak Flow Meter on.
- d) Once connected, three horizontal bars will appear on the display and last for 3 seconds. The bars will then flash to indicate that the connection between computer and device has been successfully made. As long as the cable is connected, the bars will continue to flash and the screen will remain disabled.

During the connection, the device is completely controlled by the computer. Please refer to the "Help" file in the software for detailed instructions. For more information, or help with the basic installation and use, please contact Microlife customer service at 866-934-1839

### **Notes:**

- A detailed software instruction manual is available within the software. It is accessible through the "help" button.
- Once stored memories have been downloaded from the monitor to the PC, they will be erased from the monitor's

memory.

- In order to download the data from your device, a user profile must be created (minimum required data: Family Name, Given Name, and ID).
- The ID can be any name, or combination of symbols.

## **7. Error Messages/Troubleshooting**

---

If you get a low reading this could be an early warning sign your condition is getting worse. You must follow directions exactly as instructed to obtain accurate results.

Call Microlife with questions at 866-934-1839, M-F 8 AM–8 PM Eastern Time.

If an error occurs one of the following listed error codes is displayed.

<b>Message</b>	<b>Cause</b>	<b>Remedy</b>
<b>Er1</b>	The result can not be stored in memory because date/time is not set.	Please set date and time (Section 4.2)
<b>no</b>	There is no data stored in the Monitor	Set the time and date and take a reading.
<b>Hi</b>	The result is higher than 900 ml/min.	Unit malfunction, please repeat the measurement.

## **Other possible errors and their solutions:**

If problems occur when using the device, the following points should be checked

<b>Malfunction</b>	<b>Remedy</b>
The display remains blank. Batteries have been installed.	1. Check the polarity of the batteries (+/-) 2. If the display is erratic or unusual, remove the batteries and re-install new batteries.
The instrument frequently fails to measure or the values measured are wrong.	1. Ensure that the measuring tube is connected correctly 2. Ensure that the wing wheel is rotating when air is moving through the tube. Any objects, dust liquids or mucus may interfere the rotation of the wing wheel. In this case clean the tube as explained. 3. Discuss the values with your doctor.
The unit fails to take a measurement, year is blinking	1. Ensure Date and Time are set. (See section 3.2)

If you have any questions regarding the use of this device, please ask your health care professional or contact Microlife directly at 866-934-1839. The Microlife service team will be happy to help you.


## 8. Care and Maintenance

---


### 8.1. Cleaning the Mouthpiece and the Measuring Tube

Your peak flow mouthpiece and measuring tube should be cleaned within 30 minutes after each use

1. Separate the air tube and mouthpiece from the main body. Disconnect the mouthpiece from the measuring tube by moving it in a forward direction.

2. Immerse the mouthpiece in a prepared soap solution of water and regular dishwashing soap. Thoroughly clean by swirling in the prepared soap solution. Thoroughly rinse the mouthpiece with water (distilled water is recommended). Allow the device to air dry.

**Note:** Do not tap against hard surfaces to avoid damage.



3. Now, immerse the measuring tube in a prepared soap solution and thoroughly clean by swirling in the prepared soap solution. For cleaning the measuring tube it is best to use distilled water. Use care not to damage sensitive parts. Thoroughly rinse the measuring tube with water. Allow the device to air dry.

4. Inspect the device after cleaning for signs of mucus or debris inside the air tube. Each component should be free of any film. There should be no visible debris in or on the product. The device must be visibly clean.

5. Never put the measuring tube into boiling water!

6. Reconnect the mouthpiece to the measuring tube and reassemble to the main unit as shown. When the tube is pushed completely into the fixed position, it will engage the tab on the body and click into place.

### Notes:

- If you have "hard water", mineral deposits may be visible on the mouthpiece or measuring tube. Rinse in distilled water, shake and place on a paper towel and allow to air dry before using.
- Please note temperatures lower than 50°F (10°C) and higher than 104°F (40°C) and humidity greater than 85%RH can affect the accuracy.


## 8.2. Cleaning the Main Unit

Clean the main unit once a day with a clean, damp cloth. Never put the main unit into water.

### **8.3. Maintenance and Safety**

- Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed.
- Do not use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Maintain a minimum distance of 3.3m from such devices when using this unit.
- Do not expose the device to extreme temperatures, humidity, dust or direct sunlight.
- Clean the outside of the device with a soft, dry cloth. Do not use gasoline, thinners or similar solvents. The measuring tube should not be washed in a dishwasher.
- Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.
- Never attempt to repair the instrument yourself This invalidates the manufacturer's warranty.
- If the device is dropped you may notice unusual results. Please call customer service for assistance.
- Do not use the device if you think it is damaged or if anything appears unusual.

## **9. Limited Warranty**

---

Your Digital Peak Flow Monitor is warranted for 2 years by Microlife USA Inc, against manufacturer defects for the original purchaser only, from date of purchase.

The 2 year warranty applies to the monitor. Batteries are not covered by this warranty. There are no user serviceable parts inside.

The warranty does not apply to consequential and incidental damages, or damage caused by batteries, improper handling, and accidents. Professional use, not following the operating instructions, and alterations made to the monitor or accessory by third parties, are also not included in this warranty. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

Microlife USA Inc will investigate your concern. A monitor or accessory as defined by this warranty, determined to be out of specification, will be replaced and shipped to you at no cost. A monitor or accessory as defined by this warranty, determined to be within specification, will be returned to you with a report of findings, at no cost.

Please use the customer service contact information to reach Microlife USA Inc. regarding any warranty concerns. We ask that you please contact us before sending any product back in order to better identify, and more quickly process, your concern.

## **10. Certifications**

---

<b>Device standard:</b>	IEC 60601-1 IEC 60601-1-2 IEC 60601-1-11 ATS Standard 1994 Update
-------------------------	--

### **Electromagnetic compatibility:**

Device fulfills the stipulations of the International standard IEC 60601-1-2

### **Expected Service Life:**

**Device:** 5 years

**Battery Life:** Approx. 1000 measurements

### **IP20 Classification:**

Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water. Keep Dry.

The PF100 meets the ATS recommendations for accuracy and precision in measuring peak flow.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

## **11. Technical Specifications**

---

Measuring range: .....PEF from 50 to 900 l/min.  
.....FEV1 from 0.01 to 9.99

Measuring method: .....Rotating wing wheel

Accuracy: .....ATS 1994 update  $\pm$  25 l/min or  
12% of the reading,  
whichever is greater  
FEV1  $\pm$  0.1 or  $\pm$  5% of the reading,  
whichever is greater

Measuring resolution: .....PEF: 1 l/min, FEV1: 0.01

Data safety: .....data re memorized by EEPROM

Memory: .....240 measurements with date/time

Size: .....77 (W) x 144 (L) x 48 (H) mm

Weight: .....150g (with batteries)

Storage temperature: .....-5 to +50°C or 23 to 122 °F

Storage humidity: .....10 to 90% relative humidity  
maximum

Operation temperature: .....10 to 40°C or 50 to 104 °F

Operating humidity: .....10 to 85% relative humidity  
maximum

Power source: .....2 batteries of 1.5 V, size AAA

Accessories: .....3 Mouthpieces, Storage Case

Technical alterations reserved!

Follow instructions of use to avoid hazards.

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Microlife USA, Inc.