Peak Flow Meter Replacement Mouthpiece, 3-Pack

Thank you for purchasing the Microlife mouthpiece set, suitable for use with Microlife Peak Flow meters.

Replacing the mouthpiece

- 1) Remove the old mouthpiece by pulling it away from the tube.
- 2) Clean the new mouthpiece by immersing and swirling it in a prepared solution of water and dish soap. Rinse thoroughly with distilled water and air dry.
Note: Do not tap against hard surfaces to avoid damage.
- 3) Reconnect the new mouthpiece to the measuring tube, ensuring that the curved area goes into the slot as shown.

Important: Should the measuring tube move, it must be locked back into position for the unit to work. To do this, push the tube into the main body of the meter until it clicks into place.








Care and maintenance

Your peak flow mouthpiece and measuring tube should be cleaned within 30 minutes after each use

1. Separate the air tube and mouthpiece from the main body. Disconnect the mouthpiece from the measuring tube by moving it in a forward direction.
2. Immerse the mouthpiece in a prepared soap solution of water and regular dishwashing soap. Thoroughly clean by swirling in the prepared soap solution. Thoroughly rinse the mouthpiece with water (distilled water is recommended). Allow the device to air dry.
Note: Do not tap against hard surfaces to avoid damage.
3. Now, immerse the measuring tube in a prepared soap solution and thoroughly clean by swirling in the prepared soap solution. For cleaning the measuring tube it is best to use distilled water. Use care not to damage sensitive parts. Thoroughly rinse the measuring tube with water. Allow the device to air dry
4. Inspect the device after cleaning for signs of mucus or debris inside the air tube. Each component should be free of any film. There should be no visible debris in or on the product. The device must be visibly clean.
5. Never put the measuring tube into boiling water!
6. Reconnect the mouthpiece to the measuring tube and reassemble to the main unit as shown. When the tube is pushed completely into the fixed position, it will engage the tab on the body and click into place.

Notes:

- If you have "hard water", mineral deposits may be visible on the mouthpiece or measuring tube. Rinse in distilled water, shake and place on a paper towel and allow to air dry before using.
- Please note temperatures lower than 50°F (10°C) and higher than 104°F (40°C) and humidity greater than 85% RH can affect the accuracy.

Caution

Be sure to consult the instruction booklet included with your meter for the correct use of the device. Digital instruction manuals are available on our website: www.microlifeusa.com. For use only with Microlife Peak Flow meters, not compatible with any other company's devices.

Technical specifications

Mouthpiece type:For use with Microlife Peak Flow Meters
 Weight:7.5 g per mouthpiece
 Size:30 mm diameter, 43 mm H ea
 Storage temperature:-5 to +50°C or 23 to 122 °F
 Storage humidity-10 to 90% relative humidity max.
 Operation temperature:10 to 40°C (50° to 104°F)
 Operation humidity:10 to 85% relative humidity maximum
 Expected service life:2 years (~8 months per mouthpiece)

Technical alterations reserved.
Made in China

Read the instruction manual carefully to avoid hazards.

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Microlife USA, Inc.