Premium Automatic

Arm Blood Pressure Monitor

Instruction Manual

FOR MODEL
BP3NM1-4NKRO

Important Product and Safety Information

Follow Instructions for Use. This document provides important product operation and safety information regarding this Blood Pressure Monitor. Please read this document thoroughly before using the device and keep for future reference.

This blood pressure monitor is an automatic digital blood pressure measuring device for use by adults on the upper arm at home or in the doctor's office. It enables a very fast and reliable measurement of the systolic and diastolic blood pressure, as well as the pulse rate, by way of the oscillometric method. This device detects the appearance of irregular heartbeats during measurement and provides a warning signal when the irregular heartbeat is detected.

This blood pressure monitor IS intended to be used:

- For self-measurement/monitoring of blood pressure and pulse in adults; whereas the person being measured may be the user/operator of the device.
- Within a home healthcare environment.
- With a cuff located upon the user's upper arm; ½ inch above the elbow over the artery as indicated on the cuff.

This blood pressure monitor IS NOT intended for use with:

- Pregnancy
- Preeclampsia
- Children under 12
- Neonatal patients

WARNING – Potentially hazardous situation that if not avoided may result in serious injury or death.

- Self-measuring means monitoring, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you independently alter the dosages of any

drugs prescribed by your physician.

- Consult your physician before using this device if any of the following or similar conditions are present: arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, age, pregnancy, preeclampsia, renal diseases. Motion during measurement, including trembling or shivering may affect the measurement.
- Though not for use with children under 12, ensure that any children around this device are supervised; some parts are small enough to be swallowed and any present tubes or cables may provide a risk of strangulation.
- Ensure the cuff tubing is not kinked during use as harmful injury may occur due to the effect of blood flow interference caused by high pressure in the cuff not releasing.

DO NOT:

- ✗ Use this device if you think it is damaged or anything appears unusual (e.g. sporadic operation, open packaging upon purchase).
- ✗ Use the displayed pulse for checking the frequency of heart pacemakers as this device is not suitable for this action.
- ✗ Open/modify this device; inaccuracy and/or harmful injury may result.
- ✗ Conduct many frequent measurements as this may result in harmful injury due to blood flow interference.
- ✗ Place the Cuff over a wound as this may cause additional injury.
- ✗ Place and pressurize the Cuff over/near any present intravascular access or therapy, or arteriovenous shunt, as this may cause blood flow interference and result in harmful injury.
- ✗ Place and pressurize the Cuff over a limb near the side of a mastectomy as this may cause harmful injury.
- ✗ Use the Blood Pressure Monitor on a limb simultaneously with other

medical equipment on the same limb due to possible interference with such medical equipment.

- ✗ Maintain pressure in the Cuff applied to the limb for a prolonged amount of time. Ensure that circulation in the limb is not impaired by checking circulation if prolonged/repeated exposure to pressure occurs.
- ✗ Service the device when being used or when power is supplied. When the device is not in use and power is removed, cuff assembly, batteries and AC adapter (if applicable) may be replaced by the user with Microlife supplied replacements. No other parts/components are accessible.

CAUTION – Potentially hazardous situation that if not avoided may result in minor/moderate injury, property damage, and/or damage to the device

- This device contains sensitive electronics components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g., mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring accuracy. Move the device to another location if interference is determined.
- When not using the Blood Pressure Monitor for extended lengths of time, remove the batteries to avoid potential battery leakage and damage to the monitor.
- When replacing the batteries, ensure all the batteries are replaced simultaneously to avoid battery damage and potential damage to the monitor. Microlife USA does not recommend using rechargeable batteries.
- To avoid inaccurate measurements and to lessen any discomfort from Cuff pressure, ensure the Cuff is placed correctly on the limb and fits correctly when snug (not tight), as indicated by markings with the Cuff.
- Consult your physician in cases of frequent irregular heartbeat detections.

- This Risk Indicator feature is provided in order to help you understand your potential blood pressure risk. However, this feature is neither a diagnosis nor a substitute for a medical examination. It is important to consult with your physician to determine your risk.

DO NOT:

- ✗ Drop this device or expose it to strong vibrations; sensitive components may be affected resulting in inaccuracies and/or operational issues.
- ✗ Use the Blood Pressure Monitor outside of its specified operation temperature and humidity rating, or if stored outside of its specified storage temperature and humidity rating. Avoid storage in direct sunlight.
- ✗ Use this device in a moving vehicle; inaccurate measurements may result.
- ✗ Use third party accessories. Only use Microlife authorized accessories, such as cuffs or AC adapters, as those not approved for use with the device may provide inaccurate measurements, injury, and/or damage the device.

STANDARDS

In addition to the standards stated in the Instruction Manual:

- This medical device is compliant with medical device and non-invasive blood pressure monitor standards IEC 60601-1, IEC 60601-1-2, IEC 60601-1-11, and AAMI/ANSI/IEC 80601-2-30, and
- Electromagnetic standards IEC 60601-1-2 along with FCC Part 15, and
- Clinical Testing per standard ISO 81060-2:2013 was conducted on blood pressure device using the same measurement technology.

Please note: According to international standards, your monitor should be checked for accuracy every 2 years.

TYPE

IP20 Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water. Keep Dry. 

Type BF Applied Part

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Expected Life

Monitor: 5 Years

Cuff: 2 Years

BHS

The B.H.S. (British Hypertension Society) clinical protocol was used to measure the accuracy of this product. Blood pressure units using the same measurement technology are graded “AA” for systolic/diastolic accuracy by independent investigators using the BHS protocol. This is the highest grading available for blood pressure monitors. Please see bhsoc.org for more information. (Uses the same algorithm as B.H.S. graded “AA” model number 3BT0-1)

FCC

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Changes or modifications to the product are not approved by Microlife USA and could void the user's authority to operate the equipment under FCC jurisdiction.

This equipment has been verified to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment on and off, the user is encouraged to try to correct the interference by increasing the distance between the product and the affected device; or if applicable, reorientation the receiving antenna, or moving the product's power plug to another receptacle.

Automatic Blood Pressure Monitor

Instruction Manual

Table of contents

1. Introduction

1.1. Your automatic blood pressure monitor

2. How is blood pressure measured?

2.1. What is blood pressure, and why does it fluctuate?

2.2. Measurement guidelines

2.3. What can I do to change my blood pressure?

3. Components of your blood pressure monitor

4. How do I get started?

4.1. Inserting the batteries

4.2. Setting the date and time

4.3. Tube connection

4.4. Select the measuring mode: standard or averaging

4.5. Averaging mode

5. How do I take a measurement?

5.1. Preparing to take a measurement

5.2. Common errors

5.3. Select the user

5.4. Using the cuff

5.5. Taking a measurement

5.6. Irregular Heartbeat detection

5.7. Memory: reviewing readings

5.8. Stopping a measurement

5.9. Hypertension risk indicator

5.10. Battery change indicator

5.11. Using an AC Adapter (not included)

6. Error messages/ troubleshooting

7. Care and maintenance

8. Limited warranty

9. Technical specifications

10. How to contact us

1. Introduction

1.1. Your automatic blood pressure monitor

Your blood pressure monitor is a fully automatic digital blood pressure measuring device for use by adults on the upper arm at home or in your doctor's office. It enables very fast and reliable measurement of the systolic and diastolic blood pressure as well as the pulse by way of the oscillometric method. This device offers clinically proven accuracy and has been designed to be user friendly.

Before using your blood pressure monitor, please read this instruction manual carefully and then keep it in a safe place for future reference. If you have additional questions regarding blood pressure measurements please contact your doctor.

2. How is blood pressure measured?

2.1. What is blood pressure, and why does it fluctuate?

Your level of blood pressure is determined in the circulatory center of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is affected by fine muscles in the blood vessel walls.

Your level of arterial blood pressure changes periodically during heart activity: During the "blood ejection" (Systole) the value is highest (systolic blood pressure value). At the end of the heart's "rest period" (Diastole) pressure is lowest (diastolic blood pressure value). Blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2. Measurement guidelines

Blood pressure is very high if your systolic blood pressure is over 140 mmHg and/or your diastolic pressure is above 90 mmHg. In this case, please consult your doctor immediately. Long-term values at this level endanger your health due to continual damage to the blood vessels in your body.

Should the systolic blood pressure values lie between 130 mmHg and 139 mmHg and/or the diastolic blood pressure values lie between 80 mmHg and 89 mmHg, consult your doctor. Regular self-checks will be necessary.

With blood pressure values that are too low (i.e. systolic values under 105 mmHg and/or diastolic values under 60 mmHg), consult your doctor.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your doctor. Never use the results of your measurements to independently alter the medication prescribed by your doctor.

Which values are normal?

The following standards for assessing high blood pressure in adults have been established by the American Heart Association (AHA) and American College of Cardiology (ACC) in 2017:

Category	Systolic (mm Hg)	Diastolic (mm Hg)
Normal	<120	and <80
Elevated	120-129	and <80
Hypertension		
Stage 1 Hypertension	130-139	or 80-89
Stage 2 Hypertension	140-179	or 90-119
Hypertensive Crisis	≥180	or ≥120

Additional information

- This chart reflects a 2017 update to blood pressure standards.
- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called “labile hypertension”. Consult your doctor.

2.3. What can I do to change my blood pressure?

- a) Consult your doctor.
- b) Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result from arteriosclerosis. Furthermore, the heart will become structurally damaged with increased blood pressure values.
- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organic malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take to reduce and prevent high blood pressure. These measures must be permanent lifestyle changes.

A) Eating habits

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

B) Previous illnesses

- Consistently follow all medical instructions for treating illnesses such as:
 - Diabetes (Diabetes mellitus or sugar diabetes)
 - Fat metabolism disorder
 - Gout

C) Habits

- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (coffee, tea, chocolate)




Your physical condition

- **After a medical examination, and with your doctor's approval and direction; exercise.**
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your exercise routine. Speak with your doctor regarding the type and extent of exercise that is appropriate for you.

3. Components of your blood pressure monitor

The illustration shows the blood pressure monitor consisting of:


a) Measuring unit

b) Upper arm cuffs:

For arm circumference 8.7" - 16.5"

If you ever need to buy a replacement cuff call customer service at 1-800-568-4147 for information.

Please Note:

Arm circumference should be measured with a measuring tape in the middle of the relaxed upper arm. Do not force cuff connection into the opening. Make sure the cuff connection is not pushed into the AC adapter port. If the cuff is too small, call 1-800-568-4147 for additional information. You may prefer to use one of our wrist type monitors.

4. Operation of your blood pressure monitor

4.1. Inserting the batteries

After unpacking your device, insert the batteries. The battery compartment is located on the bottom of the device.

- a) Remove cover on the bottom of the monitor.
- b) Insert the batteries (4 x size AA 1.5 V), following the indicated polarity. (+ -)
- c) If the low battery indicator appears in the display, the batteries are discharged and the unit will not function.
- d) If the batteries are inserted incorrectly, the display may function erratically or not at all. Check battery polarity.

Attention!

- After the low battery indicator appears, the device won't function until the batteries have been replaced.
- Please use AA 1.5 V batteries and replace them all at once.
- Do not use rechargeable batteries.
- If the blood pressure monitor is not used for long periods, remove the batteries from the device.


Low Battery Indicator

4.2 Setting the date and time

1. After the new batteries are activated, the year number flashes in the display. You can advance the year by pressing the "M" button. To confirm and then set the month, press the time button.
2. You can now set the month using the "M" button. To confirm and then set the day, press the time button.
3. Please follow the instructions above to set the day, hour and minutes.
4. Once you have set the last minute and pressed the time button, the date and time are set and the time is displayed.
5. If you want to change the date and time, press and hold the time button down for approx. 3 seconds until the year number starts to flash. Now you can enter the new values as described above.

4.3. Tube connection

Insert the cuff tube into the opening on the side of the monitor indicated by the drawing of a cuff.


4.4. Select the measuring mode: standard or averaging mode

This instrument enables you to select either standard (single measurement) or Averaging Mode (automatic triple measurement). To select standard mode, push the switch on the side of the instrument downward to Position "1" and to select Averaging mode, upwards to Position "3". If you select 1 then only one measurement will be taken. If you select 3, the unit will inflate and deflate three times resulting in one final average measurement.

4.5. Averaging mode

- In Averaging Mode, three measurements are automatically taken in succession and the result then automatically analyzed and displayed. Because your blood pressure constantly fluctuates, a result determined in this way is more reliable than one produced by a single measurement.
- After pressing the Start/Stop button the selected Averaging Mode appears in the display as the  symbol.
- The bottom, right hand section of the display shows a 1, 2 or 3 to indicate which of the 3 measurements is currently being taken.
- There is a break of 15 seconds between the measurements (15 seconds are adequate according to "Blood Pressure Monitoring, 2001, 6:145-147" for oscillometric instruments). A count down indicates the remaining time and a beep will sound 5 seconds before the 2nd and 3rd readings will begin.
- The individual results from measurements 1, 2 and 3 are not initially displayed.
- **If you want to see each of the individual readings making up a triplicate reading, press and hold the "M" button for 3 seconds immediately after an Averaging Mode measurement. Please note: You can only see the individual readings one time immediately following an Averaging Mode measurement. See section 5.4.**
- Do not remove the cuff between measurements.

5. How do I take a measurement?

Please note: You should always be seated before and during measurement

5.1. Preparing to take a measurement

- Avoid eating and smoking as well as all forms of exertion directly before measurement. These factors influence the measurement result. Find time to relax by sitting in an armchair in a quiet atmosphere for about 5 minutes before measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Always perform measurements at the same time of day, since the blood pressure changes during the course of the day.

5.2. Common errors:

Note: Comparable blood pressure measurements always require the same conditions. These are normally quiet conditions.

- All efforts by the user to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not flex any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower or higher than the heart, an erroneously high or low blood pressure will be measured. Each 15 cm (6") difference in height between your heart and the cuff results in a measurement error of 10 mmHg.
- Cuffs that are too small result in false measurement values. Selecting the correct cuff is extremely important. Cuff size is dependent upon the circumference of the arm (measured in the center). The permissible range is printed on the cuff. If this is not suitable for your use, please call 1-800-568-4147.
- A loose cuff or a sideways protruding air pocket causes false measurement values.

- With repeated measurements, blood accumulates in the arm, which can lead to false results. Consecutive blood pressure measurements should be repeated after a 1 minute pause to allow the accumulated blood to flow away.

5.3. Select the user

This blood pressure monitor is designed to store 99 measurements for each of two users. **Before taking a measurement, be certain that the correct user has been selected.**


- a) Press and release the TIME button. The user icon will flash.
- b) Press the MEMORY button to toggle between user 1 and user 2
- c) Press the ON/OFF button to confirm your selection.

5.4. Using the cuff


a) Pass the end of the cuff through the flat metal ring so that a loop is formed. The hook & loop material must be facing outward. (Ignore this step if the cuff has already been prepared.)

b) Slide arm through cuff until the bottom edge is about $\frac{1}{2}$ " above your elbow. Align the cuff artery mark (arrow) over the arm artery that runs down the inner side of the arm.

Important! The red strip on the edge of the cuff (Artery Mark) must lie over the artery which runs down the inner side of the arm. **TIP:** Align red artery mark to pinky finger.

c) Tighten the cuff by pulling the end and close the cuff by affixing the hook & loop material.

d) There should be little free space between the arm and the cuff. You should be able to fit 2 fingers between your arm and the cuff. Clothing must not restrict the arm. Any piece of clothing which does must be removed.

e) Lay your arm on a table (palm upward) so the cuff is at the same height as the heart. Make sure that the tube is not kinked.


f) Remain seated quietly for 5 minutes before you begin the measurement.

Comment:

If it is not possible to fit the cuff to the left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.

Cuff on right arm

5.5. Taking a measurement

After the cuff has been appropriately positioned the measurement can begin. Remain still during your measurement, do not flex muscles.

- Press the START/STOP button. The pump begins to inflate the cuff. In the display, the increasing cuff-pressure is continually displayed.
- After automatically reaching an individual pressure, the pump stops and the pressure slowly falls. The cuff-pressure is displayed during the measurement. If the inflation pressure is not sufficient, the monitor automatically re-inflates.
- When the device has detected your pulse, the heart symbol in the display begins to blink.
- When the measurement has been concluded, the cuff deflates. The measured systolic and diastolic blood pressure values as well as the pulse are now displayed.
- The measurement results are displayed until you switch the device off. The device will turn off automatically after 1 minute if no buttons are pressed.

Pumping
Pressure

Measuring

Systolic

Diastolic


Pulse

Measurement
complete

5.6. Irregular heartbeat detection

The appearance of this symbol  indicates that certain pulse irregularities were detected during the measurement. In this case, the result may deviate from your normal basal blood pressure – repeat the measurement. In most cases, this is no cause for concern. However, if the symbol appears on a regular basis (e.g., several times a week with measurements taken daily), we advise you to tell your doctor.

Please show your doctor the following explanation:

Information on frequent appearance of the irregular heartbeat symbol

This instrument is an oscillometric blood pressure monitor device that also analyzes pulse frequency during measurement. The instrument is clinically tested.

If pulse irregularities occur during the measurement, the irregular heartbeat symbol is displayed with the measurement. If the symbol appears frequently or if it suddenly appears more often than usual, we recommend the patient seek medical advice.

The instrument does not replace a cardiac examination but serves to detect pulse irregularities at an early stage.

5.7. Memory: reviewing readings

At the end of a measurement, this monitor automatically stores each result with date and time. This unit stores 99 memories for each of 2 users.

Viewing the stored values

With the unit off, press the "M" button. The display first shows "A", then shows an average of all measurements stored in the unit. Please note: Measurements for each user are averaged and stored separately. Be certain that you are viewing the measurements for the correct user.

Pressing the "M" button again displays the previous value. To view a particular stored memory, press and hold the "M" button to scroll to that stored reading.

Memory full

When the memory has stored 99 results a new measured value is stored by overwriting the oldest value.

Clear all values

If you are sure that you want to permanently remove all stored values, hold down the "M" button (the instrument must have been switched off beforehand) until "CL" appears and then release the button. If you do not want to clear the values, press the ON/OFF/START button. To permanently clear the memory, press the "M" button while "CL" is flashing.

Individual values cannot be cleared.


5.8. Stopping a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g the patient feels unwell), the START/STOP button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.

5.9. Hypertension risk indicator


The colored indicator on the left-hand edge of the display has been designed to provide a quick visual representation of your blood pressure. Once a measurement has been completed, a black bar will display onscreen next to the colored hypertension indicator. The height of the black bar will show if the measurement is within the normal (green), borderline (yellow/orange) or danger (red) range.

This classification is based on standards established by the American Heart Association (AHA) and American College of Cardiology (ACC) in 2017.

If the black bar is in the:

- **green zone**, your measurement is "Normal."
- **lower yellow zone**, it is "Elevated."
- **upper yellow zone**, it is "Stage 1 Hypertensive."
- **orange zone**, it is "Stage 2 Hypertensive."
- **lower red zone**, it is "Stage 2 Hypertensive."
- **upper red zone**, it is "Hypertensive Crisis."

5.10. Battery change indicator

Batteries almost discharged

When the batteries are approximately 75% used, the battery symbol will flash a few times as soon as the instrument is switched on (if at least one of the batteries still has some charge). Although the instrument will continue to measure reliably, you should obtain replacement batteries.

Batteries discharged – replacements required

When the batteries are discharged, the battery symbol will appear, unblinking, as soon as the instrument is switched on. You cannot take any more measurements and must replace the batteries.

1. Open the battery compartment on the bottom of the instrument.
2. Replace all of the batteries at the same time – ensure they are correctly connected, as shown on the symbols in the compartment.
3. The memory retains all values although date and time must be reset – the year number will flash automatically after the batteries are replaced.
4. To set date and time, follow the procedure described in Section 4.3.

Note:

Use four new AA 1.5 V batteries. Do not use batteries beyond their expiration date. If the monitor is not going to be used for a prolonged period, the batteries should be removed.

5.11. Using an AC adapter (not included)

You may also operate this monitor using an AC Adapter(output 6V DC/600 mA with DIN plug) - not included. Only use a Microlife compatible AC Adapter to avoid damaging the unit.

1. Ensure that the AC adapter and cable are not damaged.
2. Plug the adapter cable into the AC Adapter port on the blood pressure monitor.
3. Plug the adapter into your electrical outlet.





When the AC adapter is connected, no battery current is consumed.

Note: No power is taken from the batteries while the AC Adapter is connected to the monitor. If electrical power is interrupted, (e.g., by accidental removal of the AC adapter from the outlet) the monitor must be reset by removing the plug from the socket and reinserting the AC adapter connection.

Please consult us if you have any questions regarding AC adapters by calling 1-800-568-4147.

6. Error messages/ troubleshooting

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed.

Error No.	Possible cause(s)
ERR 1 	The systolic pressure was determined but afterwards the cuff pressure fell below 20 mmHg. This situation can occur, for example, if the tube becomes unplugged after the systolic blood pressure has been measured. Additional possible cause: No pulse has been detected.
ERR 2 	Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the measurement (Artifact).
ERR 3 	Inflation of the cuff takes too long. The cuff is not correctly seated or the hose connection is not sealed.
ERR 5 	The difference between systolic and diastolic is excessive. Measure again carefully following proper cuff procedures and ensure measurement under quiet conditions. Contact your doctor if this message occurs frequently.
HI 	The cuff pressure is too high. Relax for 5 minutes and repeat the measurement.*
LO 	The pulse is too low (less than 40). Repeat the measurement.*

***If this or any other problem occurs repeatedly, please consult your doctor.**

Other possible errors and their solutions

If problems occur when using the device, the following points should be checked:

Malfunction	Remedy
The display remains blank when the instrument is switched ON although the batteries are in place.	<ul style="list-style-type: none">• Check battery installation.• If the display is unusual, remove the batteries and then exchange them for new ones• Check polarity
The pressure does not rise although the pump is running.	<ul style="list-style-type: none">• Check the connection of the cuff tube and connect properly.
The device frequently fails to measure the blood pressure values, or the values measured are too low or too high.	<ul style="list-style-type: none">• Fit the cuff correctly on the arm. (see 5.3)• Before starting the measurement make sure that the cuff is not too tight and that clothing is not exerting pressure on the arm. Take off articles of clothing if necessary.• Measure blood pressure again in complete peace and quiet.
Every measurement results in different values, although the device functions normally and normal values displayed.	<ul style="list-style-type: none">• Please read the following information and the points listed under “Common sources of error”.(5.2)• Repeat the measurement.• Blood Pressure changes constantly. The observed readings may accurately reflect your pressure.

Malfunction**Remedy**

The blood pressure values differ from those measured by my doctor.

- Record daily values and consult your doctor.
- Pressure readings in your doctor's office may be higher due to anxiety.
- Your blood pressure fluctuates throughout the day so record the time along with your reading.

After the instrument has inflated the cuff the pressure falls very slowly or not at all. (No measurement possible).

- Check cuff connections.
- Ensure the unit has not been tampered with.

7. Care and maintenance

- a) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.
- b) Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. Do not wash the cuff in a dishwasher, clothes washer, or submerge it in water.
- c) Handle the tube carefully. Do not pull on it. Do not allow the tubing to kink and keep it away from sharp edges.
- d) Never open the monitor. This voids the manufacturer's warranty.

8. Limited warranty

Your Automatic Blood Pressure Monitor is warranted for the estimated life of the Monitor by Microlife USA Inc, against manufacturer defects for the original purchaser only, from date of purchase.

The lifetime warranty applies to the monitor only. The following accessories are warranted for 1 year: cuff. Batteries are not covered by this warranty. There are no user serviceable parts inside.

The warranty does not apply to consequential and incidental damages, or damage caused by batteries, improper handling, and accidents. Professional use, not following the operating instructions, and alterations made to the monitor or accessory by third parties, are also not included in this warranty. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

Microlife USA Inc will investigate your concern. A monitor or accessory as defined by this warranty, determined to be out of specification, will be replaced and shipped to you at no cost. A monitor or accessory as defined by this warranty, determined to be within specification, will be returned to you with a report of findings, at no cost.

Please use the customer service contact information to reach Microlife USA Inc. regarding any warranty concerns. We ask that you please contact us before sending any product back in order to better identify, and more quickly process, your concern.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

9. Technical specifications:

Operating temperature:	10 -40 °C / 50 - 104 °F
Storage temperature:	-20 to +55 °C / -4 to +131 °F
Humidity:	15 to 90% relative humidity maximum
Weight:	375 g (including batteries)
Dimensions:	89 x 123 x 57 mm
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	Oscillometric
Pressure sensor:	Capacitive
Measuring range:	
SYS:	60 to 255 mmHg
DIA:	40 to 200 mmHg
Pulse:	40 to 199 per minute
Cuff pressure display range:	0–299 mmHg
Memory:	Automatically stores 99 measurements for each user
Measuring resolution:	1 mmHg
Accuracy:	Pressure within ± 3 mmHg or 2% of reading >200 mmHg Pulse ± 5 % of the reading
Power source:	4 AA batteries, 1.5 V, AC adapter 6V DC 600 mA (not included)
Accessories:	Softfit™ cuff for arm circumference 8.7-16.5" (22-42 cm)

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