


## AC Adapter Instructions

For use with Microlife Upper Arm Blood Pressure Monitors (output 6 V DC/600 mA with DIN plug)

1. Insert the AC adapter plug into the socket in the back or side of your monitor.
2. Plug the AC adapter into a power socket. Ensure that power is connected to the monitor. If the power is interrupted during the measurement, (e.g. by accidental removal of the AC adapter from the socket) the instrument must be reset by removing the plug from its socket and re-inserting the connector.



Remember to practice electrical safety precautions while using this adapter. For additional information call the blood pressure support toll-free number: 1-800-568-4147.


**5 Year Warranty:** This AC adapter is guaranteed to be free of manufacturing defects from the date of purchase under normal use. This warranty applies to the original purchaser only. This warranty does not include any liability for incidental or consequential damages. For repairs on this AC adapter contact 1-800-568-4147.

## Instrucciones del Adaptador AC

Para utilizarse con los monitores de presión Microlife. (Salida 6 Vdc / 600 mA con conector DIN)

1. Inserte el adaptador AC dentro del conector receptor en la parte trasera ó en el lado izquierdo del monitor.
2. Enchufe el adaptador AC en un receptáculo eléctrico. Asegúrese que la potencia eléctrica está conectada al monitor. Si la potencia es interrumpida durante la medición, (ej., si se remueve el conector accidentalmente del conector receptor) el monitor debe ser reajustado. Esto se logra removiendo el adaptador del conector receptor y re-insertándolo nuevamente.

Recuerde practicar y tomar precauciones eléctricas mien-



tras utilice este adaptador. Para información adicional llame línea libre de cargo a la Líneas de Ayuda para Monitores de Presión: 1-800-568-4147.

**Garantía de por 5 años :** Este adaptador está garantizado de por vida contra defectos de manufactura. Esto aplica solamente para el comprador original, desde la fecha de compra, y bajo el uso normal del adaptador. Esta garantía no incluye ninguna responsabilidad por daños incidentales o consecuentes. Ara reparaciones de este adaptador llame al 1-800-568-4147.